
KIM System

INSTRUCCIONES DE USO
INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRAUCHSANLEITUNG
ISTRUZIONI PER L'USO

CONTASURE KIM

CE 0197

DIMA

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax : +34-976-891872
e-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
e-mail: info@neomedic.com

Date of last revision: JANUARY 2021

IMPORTANTE

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IMPORTANT

Please, read the Instructions for Use carefully.

Not following user instructions or indications may cause malfunctioning of the system or cause damage to the patient.

This document aims to explain the instructions for use of the Contasure KIM. These instructions for use are not a substitute for training to correct SUI (Stress Urinary Incontinence). Contasure KIM should only be used by doctors instructed in the treatment of this kind of pathology and more specifically in the use of this device.

These instructions for use explain the common steps for implanting the Contasure KIM. Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may be also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

INDICATIONS FOR USE

Contasure KIM is intended as a suburethral sling for the treatment of recurrent or moderate to severe female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

DESCRIPTION OF THE DEVICE AND DEVICE ELEMENTS

All components of this device have been designed for single use only. System components:

- Polypropylene monofilament mesh Sling.
- Two Passers used for transobturator mesh placement are included in system kits. The Passer needle tips are curved with a small groove to connect to the Sling. Depending on the product reference, the system may include a winged Trocar.

CONTRAINDICATIONS

Do not implant the device:

- in patients with pre-existing infections or others condition that pose an unacceptable surgical risk.
- in patients undergoing anticoagulation treatments.
- in patients who suffer urinary tract infections or obstructions.
- in pregnant women or women who are considering future pregnancies.
- in patients with renal insufficiency.
- in patients with immune diseases.
- in patients who are still growing or less than 18 years of age.
- in patients who suffers cancers of the vagina, cervix, or uterus or in patients who have received radiation in the area of treatment.
- in sub-urethral or peri-urethral areas with signs of tissue necrosis.
- in patients with known sensitivity or allergy to device materials.
- in patients undergoing concomitant bowel surgery.
- in direct contact with bowel or visceral organs, including the urinary bladder.
- in patients with other pelvic floor mesh.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to incontinence mesh surgeries, and the reason for choosing a mesh sling should be explained. Such alternatives include but are not limited to: Biofeedback, pelvic floor exercises (Kegel exercises), behavioral modification, pessary/vaginal tampons, transurethral bulking agents, medication, electrical stimulation and/or pads. The patient may also be offered non-mesh surgical options (e.g.: fascial sling, stitches).
- Conservative, non-invasive treatments should be recommended before choosing a suburethral sling implant.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implanted mesh may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The patient should be counseled that future pregnancies may negate the effects of surgical procedure and the patient may again become incontinent.
- The physician must determine the risks that should be included in the patient informed consent.
- Users should be familiar with the surgical technique for bladder neck suspensions and should be adequately trained in the placement of Contasure KIM implantation procedure before employing the device. It is important that the sling be placed without tension under the mid-urethra.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- As the anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended location for mesh placement are planned and known for each individual patient. Employment of imaging methods before and after mesh placement may aid in proper mesh placement and confirm absence of injury to non-target anatomical structures.
- The risks and benefits of performing the Contasure KIM procedure in patients with the following conditions should be carefully considered by the physician:
 - Overweight patients;
 - Concomitant pelvic organ surgery;
 - Patients with conditions that would compromise healing, and;
 - Patients with mild to light SUI.
- Contasure KIM is not recommended for women in the first two years after their last pregnancy.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomy position.
- Special care should be taken in patients with anterior prolapse because of anatomical distortion of the pelvic floor.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or reesterilize this device. No component of the Contasure KIM system is reusable. Reject and do not use any unpackaged system.
- Before surgery, remove and discard the protective tips on the Passers and, if applicable, the winged Trocar.
- The sling tension should only be adjusted intraoperatively.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of mesh infection.
- Contasure KIM mesh Sling should not be placed in contact with bowel or visceral organs, including the urinary bladder.
- Cystoscopy should be performed to confirm bladder integrity after the Passer needle passage.
- Avoid excess tensioning of the mesh to minimize urinary retention, obstruction or the possibility of damage to tissues or internal organs.
- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of needles will minimize risks.
- Do not contact any part of Contasure KIM with any staples, clips, or clamps as mechanical damage to the mesh may occur.
- No data is available on the use of Contasure KIM with tissue adhesives and is not recommended.
- The Transition Cones and the Traction Thread loops are not permanent implants.
- The discarded Transition Cones, the Traction Thread loops and excess Sling may be a potential biohazard. The Passers have sharp tips, which could potentially protrude when not in use. Handle and dispose of the Passers, as well as the discarded components in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.
- Omitting the vaginal and skin incisions to insert the Passers may result in excessive force and could cause injury to the patient.
- A digital rectal exam should be performed to detect or avoid possible rectal perforation.
- The implanted mesh sling must always be separated from the abdominal cavity by peritoneum.

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative retropubic bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- Patients must avoid sexual intercourse, more than moderate physical activities, stress or sport activities (for example, biking, jogging, etc.) for a minimum of two months following

surgery.

- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as administered to the patient for easy identification of the implant.

POTENTIAL ADVERSE REACTIONS

Known risks of surgical procedures for the treatment of urinary incontinence with a mid-urethral sling include:

- Acute inflammatory tissue reaction/swelling and transitory local irritation.
- Foreign body adverse reactions.
- Allergic reaction.
- Abscess.
- Seroma.
- Mesh erosion and extrusion/exposure in the vagina, urethra, bladder or other pelvic structures.
- Mesh contraction (shrinkage). Associate mesh contraction with vaginal shortening, vaginal tightening and vaginal pain.
- Migration of the mesh.
- Adhesion formation.
- Calcification of device.
- Infection, stimulation of a pre-existing infection.
- Fistula formation.
- Acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area. Dysuria. Other pains.
- Nerve damage, neuromuscular problems.
- Bleeding (hemorrhage, perforation of vessels, hematuria, edema, erythema, hematoma).
- Perforation of bladder, bowel, vaginal wall, urethra, ureter, intestine and other pelvic structures.
- Altered bowel function, problems defecating.
- Failure of the procedure resulting in continued or worsening urinary incontinence, recurrent urinary incontinence.
- Urgency, frequency or nocturia may appear after the intervention.
- De novo urge Incontinence may appear after the intervention.
- Excess tension on the mesh may produce temporary or permanent urinary tract obstruction and retention.
- Voiding dysfunction.
- Foreign body sensitivity or palpability, discomfort.
- Dyspareunia. All participants may experience irritation and/or pain during sexual activity.
- Atrophic vaginitis or symptoms like redness, itching, burning or dryness of the vagina.
- Vaginal spotting, bleeding or abnormal discharge.

As with all surgical intervention, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to:

- Antibiotics.
- Blood transfusions.
- Intravenous (IV) therapy.
- Drainage.
- Multiple restoration surgeries.
- Explantation of the device.

In cases in which the polypropylene mesh needs to be removed in part or whole, significant dissection may be required.

SURGICAL TECHNIQUE

Before using Contasure KIM the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical procedure to implant Contasure KIM can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

Contasure KIM placement (transobturator technique):

1.- **Vaginal Incision and Dissection:** With the patient in modified dorsal lithotomy position, prepare the surgical field with a vaginal valve. Start by making a longitudinal incision in the anterior vaginal wall, preparing the submucosal space to present the best possible view of the urethra tract. Extend the dissection laterally (periurethral) with Metzenbaum scissors. Blunt or sharp dissection may be used. Dissect the periurethral attachment to the vagina. Insert the tip of the Metzenbaum scissors laterally, spread and advance the scissors until the tip of the scissors touches the inferior portion of the inferior pubic ramus (about 1-1.5 cm). Do this bilaterally. The vaginal dissection should be large enough for a fingertip to enter in both directions. (*Note: It is suggested that the device insertion be completed on one side before beginning dissection of the second side.*)

2.- **Skin Incisions:** Before making skin incisions, palpate the edge of the ischiopubic ramus beginning at the level of the vaginal incision; continue palpating the edge of the bone moving cephalad toward the level of the clitoris and note where the adductor longus tendon inserts into the pubic ramus. Just inferior to this location, and lateral to the bony edge, will be the mark for skin incisions. Repeat on patient's contralateral side. Confirm that both marks lie in a straight line approximately at the level of the clitoris. Make small vertical incisions over these marks, over the obturator foramen.

3a.- **Inside-out:** Insertion of the Passer needle through the Retzius space: With the Traction Threads connected to the needle tip, insert the Passers from the vaginal incision through the obturator foramen and continue until the tips exit through the skin incision. The winged Trocar should be used as guidance for the consistent passage through the dissection tract, until the Passer tips appear at the skin incisions, crossing the obturator foramen through the most medial part. To do so, insert the winged Trocar into the dissected tract next to the ischio-pubic ramus and into the obturator internus muscle, but do not perforate the obturator membrane (*Image A*). The insertion zone on the winged Trocar should be visible during this part of the procedure. Placement beyond the insertion zone may allow unintended entry into the obturator membrane or space of Retzius. If bony contact is not achieved after insertion of the Trocar, remove and re-evaluate the angle of dissection. If difficulty is encountered during insertion of the Trocar, reconfirm the direction of the dissected tract with the scissors.

Insert the first Passer into the dissected tract following the channel of the winged Trocar. Ensure that the Passer handle is oriented so that the straight tip of the Passer is aligned with the channel in the Trocar and remains in that orientation until the tip reaches/contacts the superior border of the inferior ischio-pubic ramus and is pushed slightly into the internus muscle, such that the tip of the Passer is not into the obturator membrane. Once the Passer is in position (i.e., when the tip is firmly positioned behind the ischio-pubic ramus), maintain handle position while advancing the tip and radius of the Helical Passer to perforate the obturator membrane (determined by a palpable pop), then remove the winged Trocar and retain it in the sterile field for later use in the procedure. Repeat on the contralateral side.

Note: The open side of the winged Trocar must be facing the surgeon and be pressed tightly against the patient's contra-lateral thigh. The Trocar should be held in place by the surgeon or assistant.

3b.- **Outside-in:** Insertion of the Passer needle through the Retzius space: Point the needle tip perpendicular to the skin with tip in the incision. Place index finger in the vaginal incision, until fingertip palpates the bone. Rotate the needle handle, moving the needle tip and curve around the posterior surface of the ischial pubic ramus until meeting the index finger tip. Guide the Passer needle tip through the dissected tract and towards the vaginal incision. Once the tip exits the vaginal incision, insert the Traction Threads into the needle tip groove, retract Passers until the Traction Threads are pulled through to the skin incisions over the obturator foramen. Repeat on contralateral side.

4.- **Cystoscopy:** Perform a cystoscopy to determine the bladder integrity, and assure that the bladder is not perforated. If perforated, repair as necessary, withdraw the Passer and repeat step 3.


5.- **Positioning of Contasure KIM:** Place the Contasure KIM at the mid urethra with minimal to no tension, and close the vaginal incision. When ideal urethral support is achieved, cut the excess Sling protruding from the incisions, including the Traction Threads and Transition Cones. Close all incisions.

Contasure KIM removal:

The removal technique of Contasure KIM can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

- 1.- Vaginal incision: With the patient in modified dorsal lithotomy position, prepare the surgical field with a vaginal valve. Start by making a longitudinal incision in the anterior vaginal wall, preparing the submucosal space to present the best possible view of the urethral tract.
- 2.- Identify the mesh by direct view and palpation. Dissect the mesh from the surrounding structures. Avoid cutting the mesh whenever possible. Dissection of the mesh should continue by slightly pulling it to identify the mesh path. Make two groin skin incisions at the previous needle exit points of the obturator muscle. Identify the mesh exit point and dissect the mesh from the surrounding tissue. When mesh is completely dissected, remove the mesh from the vaginal side. Rinse the entire surgical site with antibiotic saline solution and close all incisions.

Connecting the Traction Thread to Passers with o-shaped groove:

MATERIALS USED IN MANUFACTURING PROCESS

Sling: Polypropylene monofilament mesh (colorless) with blue edges (solvent blue 104, CAS no.: 116-75-6)
 Sling Transition Cones: polyoxymethylene (POM)
 Sling Traction Threads: ultra-high molecular weight polyethylene (UHMWPE) multifilament
 Passers & Trocar: Stainless steel*
 Passer handle: polypropylene

No phthalates, latex or other hazardous materials are used in the manufacture of Contasure KIM mesh.
 *Stainless steel components contain nickel

STERILIZATION

Contasure KIM is sterilized with Ethylene Oxide.
 This product is a single use device and must not be re-sterilized or re-used.
 Do not use if the packaging has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light.
 Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT

Device intended to be used once only for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER		CONSULT INSTRUCTIONS FOR USE		STERILIZED USING ETHYLENE OXIDE
	BATCH CODE		USE BY DATE		STERILE BARRIER SYSTEM WITH AN ADDITIONAL PACKAGING LAYER INSIDE
	SERIAL NUMBER		MANUFACTURER		DATE OF IMPLANTATION
	DO NOT REUSE		KEEP DRY		PATIENT NAME OR PATIENT ID
	DO NOT RESTERILIZE		KEEP AWAY FROM SUNLIGHT		NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER
	DO NOT USE IF PACKAGE IS DAMAGED		DEVICE NAME		INFORMATION WEBSITE FOR PATIENTS

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