SURGICAL INSTRUMENTS

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GUIDE FOR CLEANING & STERILIZATION

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1. ANSI/AAMI ST79:2006 "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities".
2. AAMI TIR12:2010 "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."

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MANTENER SECO

IMPORTANT

Please read the Instructions For Use carefully. The surgical instruments should not be used or reused without undergoing validated cleaning/disinfection and sterilization cycles. The manufacturer recommends following standard reprocessing cycles including automated cleaning and disinfection cycles, and steam sterilization (autoclave). The surgical instruments are provided with guidelines for reprocessing (cleaning and sterilization), and include parameters which are based on the state of the art for cleaning and sterilization cycles for similar devices. Ultimately, the microbiological (sterility) and end-of-life validation are the responsibility of the end user as these guidelines have not been validated by the manufacturer.

INDICATIONS FOR USE

The reusable urogynecological surgical instruments are hand-held devices intended to aid in the surgical positioning of a pelvic floor mesh or sling manufactured by DIMA.

Note: This product is provided non-sterile. New and used instruments must be thoroughly processed according to these instructions prior to use in order to attain sterility.

WARNINGS AND PRECAUTIONS

- The surgical instrument(s) should only be used by medical specialists trained in their use. Surgeons should not begin clinical use of the surgical instrument without reviewing the Instructions for Use. It is imperative that appropriate personnel read all guideline materials prior to surgery. If the guidelines are not properly followed, the risk of complications to the patient increases.
- Universal precautions should be observed by all hospital personnel that might be working with the contaminated or potentially contaminated surgical devices. Personal protective equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated medical devices and equipment. PPE should include, among other options: gown, facial mask, goggles or face shield, gloves and shoe covers in accordance with applicable health and safety procedures.
- Metal brushes or scouring pads must not be used during manual cleaning procedures. These materials will damage the surface and finish of instruments. Soft-bristled, nylon brushes and pipe cleaners should be used.
- Cleaning agents with low foaming surfactants should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument below the surface of the cleaning solution to prevent formation of aerosols and splashing which may spread contaminants. Cleaning agents must be easily and completely rinsed from device surfaces to prevent accumulation of detergent residue.
- Do not allow contaminated devices to dry prior to reprocessing. All cleaning and sterilization steps are facilitated by not allowing blood, body fluid, tissue debris, saline, or disinfectants to dry on used instruments.
- The instrument must be manually cleaned immediately after use. This initial cleaning should not be postponed until the final cleaning and resterilization procedure.
- Saline and cleaning/disinfecting agents containing aldehyde, mercury, active chlorine, bromine, bromide, iodine or iodide are corrosive and should not be used; instruments must not be placed or soaked in Ringers Solution
- Descaling agents that include morpholine should not be used in steam sterilizers. These agents leave residue which can damage polymer instruments over time.
- As with any surgical instrument, careful attention should be made to assure that excessive force is not placed on this instrument. Excessive force can result in failure.
- If the instrument does not function properly, immediately discontinue use.
- Follow instructions and warnings as issued by manufacturers of any decontaminants, disinfectants and cleaning agents used. Wherever possible, avoid the use of mineral acids and harsh, abrasive agents.
- No part of the process shall exceed 145°C. Do not sterilize by dry heat.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions for Use.

MATERIALS USED IN THE MANUFACTURE OF THE INSTRUMENT(S)

- Stainless steel
- Polypropylene - *in the case that the instrument(s) contain non-metallic parts*
- Polyphenylsulfone (PPSU) - *only the Anchoring Handle*

LIMITATIONS AND RESTRICTIONS ON REPROCESSING

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

If instruments are determined to reach end of life, they must be properly disposed of. Any used components without reprocessing may be a potential biohazard. Handle and dispose of any surgical instruments in accordance with any local and federal law regarding medical waste, with special attention to disposal requirement for sharp devices, to prevent any unnecessary injury or harm.

CARE AND MAINTENANCE INSTRUCTIONS

Proper care and maintenance is important for the efficient and safe operation of medical surgical equipment. We recommend careful inspection of all equipment upon receipt and prior to each use, as a safeguard against possible injury to patient and/or user (operator).

- Visually inspect the surgical instrument(s) prior to and after cleaning for damage or wear. Discard surgical instrument(s) with damage or excessive wear.
- The tips and groove tips (if applicable) should be smooth without denting, cracks, burrs or fractured surfaces.
- All other surfaces should be uniform and smooth, without corrosion and dents, sharp edges or points. Minor surface hazing and scratches on the polymer handles or components are acceptable. The metallic tip or needle should be firmly attached to the polymer handle or component.
- Labeling should be legible.
- If any defects are observed, contact the manufacturer. Refer to contact information contained in these Instructions for Use.

PREPARATION FOR USE

First use

Before first use, the instrument(s) must be cleaned and sterilized.

Both a manual and an automatic cleaning method are proposed in these instructions for use.

A sterilization process is proposed in these instructions for use.

Subsequent uses

Surgical instrument(s) can be reused after conducting appropriate cleaning and sterilization.

Both a manual and an automatic cleaning method are proposed in these instructions for use.

A sterilization process is proposed in these instructions for use.

Important note: The instrument must be manually pre-cleaned immediately after use. A pre-cleaning method is proposed in these instructions for use.

CLEANING INSTRUCTIONS

Both manual and automatic cleaning require a pre-cleaning.

Pre-cleaning procedure

To prevent the drying of blood and tissue residues and growth of microorganisms, instruments should be immersed, as soon as possible after each procedure, in an aqueous solution and preferably a combined neutral or enzymatic cleaning and/or disinfectant solution prepared and maintained per instructions of the solution manufacturer. This immersion process should be carried out as soon as reasonably practical following use of instrument(s).

Regardless, the cleaning and re-sterilization process of the instrument(s) should also be performed as soon as reasonably possible.

*Note: Surgical instruments with detachable components **must be** disassembled before undergoing pre-cleaning and cleaning cycles.*

Manual cleaning procedure

Upon termination of use, thoroughly rinse instrument(s) under cool running tap water (~20°C/68°F) for at least 90 seconds to remove all organic residue and contaminants. Use of a non-metallic soft bristle brush for light scrubbing is acceptable, paying special attention to areas where contaminants may be trapped (corners, grooves and holes.) A syringe to flush the lumen and hard to reach areas of the surgical instrument(s) is recommended.

- For 10 minutes, fully immerse the instrument(s) in enzymatic detergent solution for stainless steel (following the instructions for use of the detergent manufacturer) using lukewarm (approximately ~33°C/91°F) tap water. Verify that the expiration date of the detergent has not passed.
- Thoroughly clean the instrument, scrubbing with a non-metallic bristle brush. This part of the cleaning process should be executed with the instrument(s) fully submerged in the enzymatic solution to avoid splashing, paying special attention to areas where contaminants may be trapped (corners, grooves, crevices and/or lumens).
- Rinse the instrument(s) thoroughly in cold (approximately 20°C/68°F), running tap water.
- For at least 10 minutes, fully immerse the instrument(s) in a disinfectant detergent solution (following the instructions for use of the detergent manufacturer) and lukewarm (approximately 33°C/91°F) running tap water.
- Thoroughly clean the instrument, scrubbing with a non-metallic bristle brush. This part of the cleaning process should be executed with the instrument(s) fully submerged in the detergent solution to avoid splashing, paying special attention to areas where contaminants may be trapped (corners, grooves, crevices and/or lumens). A syringe to flush the lumen and hard to reach areas of the surgical instrument(s) is recommended.
- Rinse the instrument(s) thoroughly with an abundant amount of distilled water, paying special attention to the hard-to-reach areas where contaminants may remain (corners, grooves, crevices and/or lumens).
- Visually inspect the instrument for cleanliness, paying special attention to the hard-to-reach areas where contaminants may remain (corners, grooves, crevices and/or lumens).
If any residue is detected, return to STEP 1 and repeat the cleaning process completely. Once repeated, if no residue is detected, continue on with the drying procedure.
- Instrument(s) must be dried as soon as possible after cleaning in a protected or clean environment to minimize new contamination.
For automatic drying cycles, dry at 93°C (199°F) for 30 minutes.
For manual drying, the use of a swab or pressurized air in and on all surfaces and channels to dry is recommended. A single-use, lint-free, non-woven, soft cloth may be used to dry instruments.
- After drying, inspect for condensation and if any is found, repeat the drying process until no condensation is observed, then proceed with sterilization.

Automatic cleaning procedure

Note: A low-foaming, PH neutral detergent is recommended for automated cleaning cycles.

Instrument(s) may be cleaned using ultrasounds. Hereafter a possible cleaning cycle is described:

- Insert the instrument(s) into the ultrasonic washing machine.
- Fill the tank of the ultrasonic machine with detergent (following the instructions for use of the detergent manufacturer) and running tap water with a temperature of 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
- Fill the tank of the ultrasonic machine with fresh detergent (following the instructions for use of the detergent manufacturer) and running tap water with a temperature of 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
- Fill the tank of the ultrasonic machine with distilled water at a temperature between 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank.
- Fill the tank of the ultrasonic machine with fresh distilled water at a temperature between 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
- Visually inspect the instrument for cleanliness, paying special attention to the hard-to-reach areas where contaminants may remain (corners, grooves, crevices and/or lumens).
If any residue is detected, return to STEP 1 and repeat the cleaning process completely. Once repeated, if no residue is detected, continue on with the drying procedure.
- Instrument(s) must be dried as soon as possible after cleaning in a protected or clean environment to minimize new contamination.
For automatic drying cycles, dry at 93°C (199°F) for 30 minutes.
For manual drying, the use of a swab or pressurized air in and on all surfaces and channels to dry is recommended. A single-use, lint-free, non-woven, soft cloth may be used to dry instruments.
- After drying, inspect for condensation and if any is found, repeat the drying process until no condensation is observed, then proceed with sterilization.

STERILIZATION PROCEDURE

*Note: Surgical instruments with detachable components **must be** disassembled before undergoing sterilization cycles.*

- Instrument(s) will be properly packaged in order to guarantee the durability of the sterilization procedure.
- Once instrument(s) are correctly packaged, they will be inserted into an autoclave to proceed with steam sterilization.
- The person responsible for carrying out the sterilization process should carefully follow the instructions for use provided by the autoclave manufacturer, given that adequate times and temperatures to guarantee complete sterilization of instrument(s) can vary and depend on the indications outlined by each specific manufacturer.
- Once the autoclave sterilization process has finished, the packaged instrument(s) will be removed from the autoclave and will need to cool down.
- It is absolutely necessary to verify that the packaging of the sterilized instrument(s) is completely dry. If condensation is observed in or on the packaging, they are considered unsterile and must be reprocessed.
- Once the instrument(s) cool, they will then proceed to be stored in a place that guarantees the complete preservation of the packing and the instrument(s).

The operator should follow the instructions for steam sterilization provided by the sterilizer's manufacturer which indicate cycle parameters. The following sterilization times and temperatures are recommended for pre-vacuum autoclave steam sterilization:

TEMPERATURE	TIME
132°C/270°F	4-5 minutes
135°C/275°F	3 minutes

Drying time: 30 minutes

The previously described procedure should be considered as a general recommendation given that the sterilization procedure should be meticulously established following the instructions provided by the autoclave manufacturer.

STORAGE

Store wrapped surgical instrument in a designated sterile supply area that protects sterile items and their packaging from damage. Recommended storage conditions: in a dry place at room temperature, away from heat and direct light.

RESPONSIBILITIES

DIMA SL as manufacturer and distributor of the instruments, which are the subject of this manual, is not responsible for any direct or indirect damages resulting from improper or incorrect use, care, or servicing of these instruments, or from the failure to follow the guidelines contained in this manual.

If non-authorized persons perform any repairs on or modifications of this product, DIMA SL does not accept any liability and the warranty becomes void.

This product and manual are subject to technical changes and modifications without notice.

It is recommended to follow the instructions for use provided by the manufacturer. It is the responsibility of the processor to ensure that the reprocessing is performed using proper equipment, materials and personnel to achieve the desired result. Any deviations from these instructions by the processor should be validated prior to use for effectiveness and potential adverse consequences. Ultimately, the end-user is responsible for the reprocessing and safe re-use of the device, including microbiological validation.

The responsibility of developing a proper sterilization procedure for the instrument(s) lies with the organization which establishes the procedure, *not* with the manufacturer of the instrument(s).

ADDITIONAL INFORMATION

Precautions with respect to the transmission of Creutzfeldt-Jacob disease (CJD)

- CJD is the human transmissible spongiform encephalopathy (TSE) variant agent which is recognized as resistant to the normal disinfection and sterilization procedures.
- The intended use for instrument(s) is recognized as low risk procedures for the transmission of CJD. Deactivation of the CJD agent is being actively researched at this time and is beyond the scope of these instructions
- Facilities reprocessing surgical instruments must take particular precautions when handling instruments used on patients at high risk for CJD to not only deactivate the CJD agent but also to reduce cross-contamination during the transportation and cleaning processes.

References

These recommendations are developed using the following guides:

1. ANSI/AAMI ST79:2006 "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities".

2. AAMI TIR12:2010 "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."

Note: All cleaning and sterilization procedures should be validated by the health care facility performing these operations. The recommendations included in this cleaning and sterilization guide for reusable surgical instruments manufactured by DIMA are supported by the following entities:

1. Department of health and social care – GOV.UK in the following guideline: "Health Technical Memorandum 01-01 on the management and decontamination of surgical instruments (medical devices) used in acute care" published March 20, 2013 / updated July 8, 2016: www.gov.uk/government/publications/management-and-decontamination-of-surgical-instruments-used-in-acute-care
2. U.S. Department of Health and Human Services, Food and Drug Administration (FDA) in the following guideline: "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling" published March 17, 2015 / updated June 9, 2017: <https://www.fda.gov/media/80265/download>

SYMBOLS USED ON LABELING

NON-STERILE

DO NOT USE IF PACKAGE IS DAMAGED AND
CONSULT INSTRUCTIONS FOR USE

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1. Department of health and social care – GOV.UK nelle seguenti linee guida: "Health Technical Memorandum 01-01 on the management and decontamination of surgical instruments (medical devices) used in acute care" pubblicate il 20 marzo 2013/aggiornate l'8 luglio 2016: www.gov.uk/government/publications/management-and-decontamination-of-surgical-instruments-used-in-acute-care
2. U.S. Department of Health and Human Services, Food and Drug Administration (FDA) nelle seguenti linee guida: "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling" pubblicate il 17 marzo 2015/aggiornate il 9 giugno 2017: <https://www.fda.gov/media/80265/download>

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