REMEEEX READJUSTMENT KIT

REMEEEX READJUSTMENT KIT LONG

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
GEBRAUCHSANWEISUNG
MODE D'EMPLOI
ISTRUZIONI PER L'USO
INSTRUÇÕES DE USO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

NEOMEDIC INTERNATIONAL SL
C/ Maestrat 41-43, 08225 Terrassa (Barcelona), Spain
+34-93-7804505, info@neomedic.com, www.neomedic.com

DESARROLLO E INVESTIGACIÓN MÉDICA ARAGONESA S.L. (DIMA SL)
Polígono Industrial Mediavega, parcela 2.8, 50300 Calatayud (Zaragoza), Spain
+34-976-884065, dimasl@dimasl.com, www.dimasl.com

REMEEX READJUSTMENT KIT
REMEEX READJUSTMENT KIT LONG

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
GEBRAUCHSANWEISUNG
MODE D'EMPLOI
ISTRUZIONI PER L'USO
INSTRUÇÕES DE USO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

IMPORTANT

Please, read the Instructions for Use carefully.

Not following user instructions or indications may cause malfunctioning of the system or cause damage to the patient.

These instructions for use explain the common steps for use of a Remeex Readjustment kit. These instructions for use are not a substitute for training to correct SUI (Stress Urinary Incontinence). The device should only be used by doctors after undertaking specialized training, and with extensive experience in the treatment of pelvic floor disorders, more specifically in the use of this device, and by following the instructions for use.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective use of the device. As with any surgical intervention, special circumstances may also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

The Eudamed public website can be accessed at the following address: <https://ec.europa.eu/tools/eudamed>. Information for Remeex Readjustment kits may be found by performing a search of its Basic UDI-DI: 843656835REMEEXREADJUSTEE.

INDICATIONS FOR USE

Remeex Readjustment Kit is a mechanical system that permits readjustment of the Remeex System sling tension level and urethral support throughout the patient's lifetime.

DESCRIPTION OF THE DEVICE

This device is provided sterile, and has been designed for single-use.

All Remeex Readjustment Kits are accessories for Remeex System devices. The Remeex Readjustment kits are composed of 3 readjustment tools for post-operative regulation of a Remeex SUI Sling: one Manipulator and two Disconnectors. These instructions for use are valid for the following commercial references of the Remeex Readjustment kits: AR-01; AR-01LP.

CONTRAINDICATIONS

Only use the device on patients who have received a Remeex System implant.

Do not use the device:

- in patients with pre-existing conditions that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- if the fascia of the rectus abdominal muscles is damaged or does not offer enough resistance.
- in patients with known sensitivity or allergies to product materials.
- if the wound incision from the previous readjustment has not healed completely.

WARNINGS AND PRECAUTIONS

Preoperative:

- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications.
- Patient counseling should include a discussion that the device is used in conjunction with a permanent implant, and that some complications associated with the implant may require additional surgery that may or may not correct the complication.
- The physician must determine the risks that should be included in the patient informed consent.
- Users must be familiar with the surgical technique for suburethral slings and must be adequately trained in the readjustment procedure for the implanted Remeex System before employing the device.
- Surgical risks not associated directly with the readjustment, require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or resterilize this device. No component of the Remeex Readjustment kits is reusable. Reuse of the device could lead to potentially serious consequences for the patient that may lead to medical intervention, permanent injury or death. Reject and do not use any unpackaged system.
- Sling tension-level readjustments may be made throughout the lifetime of the patient. However, wound dehiscence and/or infection may occur if the surgeon performs an incision for readjustment before the previous incision wound has healed completely.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product requires a minor surgical procedure; strict aseptic measures must be taken.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of product infection.
- Avoid excess tension of the mesh to minimize urinary retention, bowel obstruction or the possibility of damage to tissues or internal organs.
- Intra-operative bleeding may occur.
- Take special care during the surgical intervention to avoid damage of nerves, blood vessels, or mechanical disruption of the tissue.
- Take special care during the surgical intervention to avoid damage to the Traction Threads which connect the Sling to the Varitensor. If the Traction Threads are damaged, treatment failure is probable.
- The use of electrocautery is not recommended during the readjustment procedure given that the device material is sensitive to excessive temperatures and may cause further risks to the patient, product and/or healthcare specialist.
- Do not turn the Manipulator without complete control. It is crucial to record the number of turns, the direction of the turns (both clockwise and counter-clockwise) and the date that it has been performed.
- While the external Manipulator is connected to the Varitensor, take the appropriate aseptic techniques and conditions to avoid contamination and infection at the site. The appropriate cleaning of the abdominal site must be performed by the appropriate medical personnel. During this time when the Manipulator and Varitensor are connected, the abdominal area must be protected using aseptic technique while covering the surgical field with appropriate sterile materials.
- To disconnect the external Manipulator from the Varitensor, gently pull on the Manipulator while rotating the Disconnector until the point mark of the Disconnector is aligned with the line mark of the Manipulator. Do not apply excessive pressure during this procedure.
- Manipulator and Disconnectors are not permanent implants.
- The Manipulator must not be connected to the Varitensor more than 3 consecutive days.
- The discarded device components (Manipulator and Disconnectors) may be a potential biohazard. Handle and dispose of the discarded components in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths, recurrent incontinence, voiding dysfunction, problems defecating or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Traction Threads may fracture in case of overload due to overexertion of the patient.

- Sling tension-level readjustments may be made throughout the lifetime of the patient. However, wound dehiscence and/or infection may occur if the surgeon performs an incision for readjustment before the previous incision wound has healed completely.
- The Physician(s) and patient should report any serious incident in relation to Remeex System to the manufacturer and the national health authority having jurisdiction in their locale.

POTENTIAL ADVERSE REACTIONS

- Acute inflammatory tissue reaction/swelling and transitory local irritation, allergic reaction or other adverse reactions to foreign bodies may occur.
- Infection, inflammation/swelling or edema, recurrence of the incontinence, discomfort and/or pain, dysuria, bleeding, erythema, hematomas, hemorrhage, abscess, and/or seromas.
- As a foreign body, the device could stimulate a pre-existing infection.
- Failure of the procedure resulting in continued or worsening urinary incontinence.
- Voiding dysfunction/de novo incontinence in the form of urgency, frequency or nocturia may appear after the readjustment.
- Excess tension on the mesh may produce temporary or permanent urinary tract obstruction and retention.
- Altered bladder and/or bowel function. Problems defecating.
- Acute and/or chronic pain.
- Nerve damage, neuromuscular problems.
- Dyspareunia. Sexual partners may also experience discomfort, irritation and/or pain during sexual activity.

As with all surgical intervention, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

In cases in which the Varitensor, and/or the Sling need to be removed in part or whole, significant dissection may be required.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, removal of the device.

SURGICAL TECHNIQUE

Before using the system, the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical procedure to readjust the Remeex System can be performed under local anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

Surgical technique

Medium-Long term ability to readjust the sling support level:


(For medium-long term readjustment, wait at least 30 days from hospital discharge to ensure that the patient is completely recovered from surgery).

Remeex Readjustment Kits offer the possibility to modify the Sling support level whenever needed throughout the patient's lifetime. During the month following the surgical procedure the tissue inflammation will diminish considerably, the manipulated organs will return to their normal anatomical position and the patient will return to her normal physical activity level. Some patients may show some degree of incontinence, in which case the Sling can be readjusted to restore continence by performing a minor surgical procedure. This is accomplished using the single-use, pre-sterilized Remeex Readjustment Kit which contains 1 Manipulator and 2 Disconnectors.

Sling Level Readjustment: In an appropriate surgical environment, after injecting local anaesthesia, make a minimal incision at the abdominal midline, just over the pubis in order to localize the Varitensor. Connect the sterile Manipulator to the implanted Varitensor and close the incision with the Manipulator protruding from the abdomen. The Sling readjustment is performed with the patient in standing position, performing the normal pressure maneuvers that cause the patient to become incontinent. Fill the bladder with 300 c.c. of saline through a urethral catheter and ask the patient to stand up and perform Valsalva maneuvers. If necessary, rotate the Manipulator clockwise checking continence level every 4 complete turns, until continence is achieved (there is no leakage). Invite the patient to urinate and measure the residual in the bladder after urination. If residual is less than 100 c.c. disconnect the Manipulator from the Varitensor, close the incision and discharge the patient. If residual is more than 100 c.c., decrease the Sling tension by rotating the Manipulator counterclockwise, helping to descend the urethra using a rigid probe. To disconnect the Manipulator from the Varitensor, insert the Disconnector inside the Manipulator and rotate the Disconnector $\frac{1}{4}$ of a turn in relation to the Manipulator while slightly pulling the Manipulator. Close the abdominal incision.

MECHANICAL FUNCTION STEPS:

- Figure 1.- To attach the Manipulator to the Varitensor, insert the Disconnector into the Manipulator, then turn it 90° in relation to the Manipulator. Remove the Disconnector (The Manipulator should be joined to the Varitensor).
- Figure 2.- Rotate the Manipulator to increase or decrease Sling tension level.
- Figure 3.- To detach the Manipulator from the Varitensor, insert the Disconnector into the Manipulator, then turn it 90° in relation to the Manipulator. Remove the Disconnector and separate the Manipulator from the Varitensor before completely closing the abdominal incision.

MECHANICAL CHARACTERISTICS /SPECIFICATIONS

Manipulator function: Every 360° turn of the Manipulator causes an increase or decrease of the Sling's suspension level by approximately 1mm.

Wind up capacity: The wind up capacity of the Varitensor is 15 cm (Female Remeex System) of Traction Threads per each side.

Maximum traction stress: The Varitensor can support the force of 18 Newtons produced by the Traction Threads, with a minimum safety coefficient of 5.

MATERIALS USED IN MANUFACTURING OF MANIPULATORPatient-contacting materials

Product references AR-01, AR-01LP: Ultra high molecular weight polyethylene (UHMWPE) + titanium

No phthalates, latex or other hazardous materials are used in the manufacture of The Remeex Readjustment kit components.

STERILIZATION

Remeex Readjustment Kits are sterilized with Ethylene Oxide. This product is a single use device and must not be re-sterilized or re-used. Do not use if the package has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and sunlight. Do not use after the expiration date detailed on the package.

SINGLE-USE PRODUCT DISCLAIMER

Device intended to be used only once and for a single patient. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER		DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		KEEP DRY
	SERIAL NUMBER		DISTRIBUTOR		DATE OF MANUFACTURE
	UNIQUE DEVICE IDENTIFIER		CONSULT INSTRUCTIONS FOR USE		KEEP AWAY FROM SUNLIGHT
	MEDICAL DEVICE		USE BY DATE		STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE		MANUFACTURER		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE
	DO NOT RESTERILIZE				

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