
CONTASURE NEEDLELESS

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MODE D'EMPLOI
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO

neomedic
international

NEOMEDIC INTERNATIONAL SL
C/ Maestrat 41-43, 08225 Terrassa (Barcelona), Spain
+34-93-7804505, info@neomedic.com, www.neomedic.com

DIMA

DESARROLLO E INVESTIGACIÓN MÉDICA ARAGONESA S.L. (DIMA SL)
Polígono Industrial Mediavega, parcela 2.8, 50300 Calatayud (Zaragoza), Spain
+34-976-884065, dimasl@dimasl.com, www.dimasl.com

CE 0197

IMPORTANT

Please, read the Instructions for Use carefully.
Not following user instructions or indications may cause malfunctioning of the system or cause damage to the patient.

These instructions for use explain the common steps for implanting Contasure Needleless. These instructions for use are not a substitute for training to correct SUI (Stress Urinary Incontinence). The device should only be used by doctors instructed in the treatment of this condition, more specifically in the use of this device, and by following the instructions for use.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may be also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

The device summary of safety and clinical performance (SSCP) is made publically available through the European Database on Medical Devices (Eudamed), and provides healthcare workers and relevant patients critical information regarding the medical device in order to enhance the reader's understanding of the device's clinical safety and performance. The Eudamed public website can be accessed at the following address: <https://ec.europa.eu/tools/eudamed>. The SSCP for the device Contasure Needleless may be found by performing a search of its Basic UDI-DI: 843656835NEEDLELESSJJ.

INDICATIONS FOR USE

Contasure Needleless is indicated as a suburethral sling for the treatment of recurrent or moderate to severe female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

DESCRIPTION OF THE DEVICE AND DEVICE ELEMENTS

All components of this device have been designed for single use only. Contasure Needleless system components:

- Polypropylene monofilament mesh Sling.
- Insertion Tip for guided passage of the Sling during placement.

CONTRAINDICATIONS

Do not implant the device:

- in patients with pre-existing infections or others condition that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- in patients undergoing anticoagulation treatments.
- in patients who suffer urinary tract infections or obstructions.
- in pregnant women or women who are considering future pregnancies.
- in patients with renal insufficiency.
- in patients with immune diseases.
- in patients who are still growing or less than 18 years of age.
- in sub-urethral or peri-urethral areas with signs of tissue necrosis.
- in patients with known sensitivity or allergy to device materials.
- in patients undergoing concomitant bowel surgery.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to incontinence mesh surgeries, and the reason for choosing a mesh sling should be explained. Such alternatives include but are not limited to: Biofeedback, pelvic floor exercises (Kegel exercises), behavioral modification, pessary/vaginal tampons, transurethral bulking agents, medication, electrical stimulation and/or pads. The patient may also be offered non-mesh surgical options (e.g.: fascial sling, stitches).
- Conservative, non-invasive treatments should be recommended before choosing a suburethral sling implant.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implanted mesh may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The patient should be counseled that future pregnancies may negate the effects of surgical procedure and the patient may again become incontinent.
- The physician must determine the risks that should be included in the patient informed consent.
- Users should be familiar with the surgical technique for bladder neck suspensions and should be adequately trained in the placement of Contasure Needleless implantation procedure before employing the device. It is important that the sling be placed without tension under the mid-urethra.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- As the anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended location for mesh placement are planned and known for each individual patient. Employment of imaging methods before and after mesh placement may aid in proper mesh placement and confirm absence of injury to non-target anatomical structures.
- The risks and benefits of performing the Contasure Needleless procedure in patients with the following conditions should be carefully considered by the physician:
 - Patients with other pelvic floor meshes;
 - Overweight patients;
 - Concomitant pelvic floor surgery, including additional SUI treatment(s);
 - Patients with conditions that would compromise healing (in patients who suffer cancers of the vagina, cervix, or uterus or in patients who have received radiation in the area of treatment, etc.) and;
 - Patients with mild to light SUI.
- Contasure Needleless is not recommended for women in the first two years after their last pregnancy.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomy position.
- Special care should be taken in patients with anterior prolapse because of anatomical distortion of the pelvic floor.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or resterilize this device. No component of the Contasure Needleless system is reusable. Reuse of the device could lead to potentially serious consequences for the patient that may lead to medical intervention, permanent injury or death. Reject and do not use any unpackaged system.
- The sling tension should only be adjusted intraoperatively.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of mesh infection.
- Avoid excess tension of the mesh to minimize urinary retention/obstruction, bowel obstruction or the possibility of damage to tissues or internal organs. Mesh must not be placed in direct contact with bowel or visceral organs, including urinary bladder.
- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper use of surgical instruments will minimize risks.

- Do not contact any part of Contasure Needleless sling with any staples, clips, or clamps as mechanical damage to the mesh may occur.
- No data is available on the use of Contasure Needleless with tissue adhesives and is not recommended.
- Omitting the vaginal incisions to insert the Insertion Tip and Sling may result in excessive force and could cause injury to the patient.
- The Sling central traction thread and the Insertion Tip are not permanent implants.
- If used, be sure to withdraw the Insertion Tip.
- Cut and withdraw the blue, central thread before closing the vaginal wall incision.
- The discarded Insertion Tip and Sling central Traction Thread may be a potential biohazard. The Insertion Tip could potentially protrude when not in use. Handle and dispose of the Insertion Tip, as well as the discarded components in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- It is recommended that patients do not overly exert themselves physically (e.g.: heavy lifting) or partake in any moderate to high intensity sporting activity (e.g.: biking, jogging, etc.) for at least a minimum of two months. It is also recommended that patients avoid sexual intercourse during the first two months after surgery.
- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as an implant card to be administered to the patient for easy identification of the implant.
- Taking into account the condition of each individual patient, the surgeon will plan for post-operative follow-up. Follow-up is recommended at short/mid-term and long-term if possible.
- The Physician(s) and patient should report any serious incident in relation to Contasure Needleless to the manufacturer and the national health authority having jurisdiction in their locale.

POTENTIAL ADVERSE REACTIONS

Known risks of surgical procedures for the treatment of urinary incontinence with a mid-urethral sling include:

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur.
- Inflammation/swelling, adhesion formation, fistulas, abscess, seromas and/or hardening of the mesh caused by fibrosis.
- Bleeding (hemorrhage, perforation of vessels, hematuria, edema, erythema, hematoma).
- Perforation of bladder, bowel, vaginal wall, urethra, ureter, intestine and other pelvic structures.
- As a foreign body, the device could stimulate a pre-existing infection, or suffer rejection.
- Extrusion, exposure, and/or migration of the mesh from the desired location. The device may also produce erosion through the vagina, urethra, bladder or other surrounding tissues or structures. Mesh contraction (shrinkage) is associated with vaginal tightening, vaginal shortening and vaginal pain.
- Urinary problems such as voiding symptoms/de novo incontinence in the form of urgency/overactive bladder (OAB), frequency or nocturia may appear after the intervention.
- Failure of the procedure resulting in recurrence (continued or worsening urinary incontinence).
- Altered bladder and/or bowel function such as temporary or permanent urinary retention/obstruction or problems defecating.
- Infection, urinary infection.
- Acute and/or chronic pain in the groin, thigh, buttox, leg, pelvic and/or abdominal area. Dysuria. Other pains. Discomfort.
- Nerve damage, neuromuscular problems.
- Vaginal spotting, bleeding, abnormal discharge or scarring.
- Dyspareunia. Sexual partners may experience discomfort, irritation and/or pain during sexual activity.
- Foreign body sensitivity or palpability.

As with all surgical intervention, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

In cases in which the mesh needs to be removed in part or whole, significant dissection may be required.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries and/or removal of the device.

SURGICAL TECHNIQUE

Before using Contasure Needleless the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical procedure to implant Contasure Needleless can be performed under local, regional or general anesthesia. It is recommended to follow the surgical implants antibiotic protocol dictated by the hospital.

Contasure Needleless placement:

With the patient in lithotomic position, inject local anaesthesia at the anterior vaginal wall.

Make a 2 cm longitudinal incision at the middle line of the anterior vaginal wall. Variations may occur in specific incisions due to individual technique and patient anatomy.

Dissect the para-urethral submucosa with Mayo scissors until the tip of the scissors is in contact with the inferior pubic ramus. Do not dissect the fascia with Mayo scissors. Expand the submucosal space through wide retrodissection while withdrawing the Mayo scissors to create enough space for the mesh to be absolutely flat under the urethra. Insert surgical forceps (Crawford) into one of the terminal pockets of the Sling, open the forceps inside the pocket to gather mesh, closing them to fold the pocket, and close it to minimize the width of the sling pocket*. Insert the sling through the submucosal dissection in the direction of 10 o'clock for the right side of the patient and 2 o'clock for the left side of the patient. Push the mesh pocket through the fascia into the internal obturator muscle until to the point where the central sling traction thread is 0.5 to 1cm from the urethral midline. Open the surgical forceps to release the mesh pocket. Close and pull the forceps out of the vaginal surgical incision maintaining the mesh in place. Repeat the same maneuver on the contra lateral side and continue with forceps penetration up to the point where the central traction thread of the mesh is at the midurethral line. A Mosquito clamp can be placed on the central traction thread to better control the center of the sling.

*Note: Be certain that the entire Sling stays within the sterile field during passage through the tissues.

Sling repositioning maneuvers can be performed by further pushing the sling pocket with the tip of the forceps. Ensure that the sling is fully extended, in full contact with the submucosa and centered at the mid urethral line.


Cut and withdraw the traction thread at the center of the sling. If Insertion Tip has been used during Sling placement, be sure to completely withdraw it before closing the vaginal incision. Close the vaginal incision.

Contasure Needleless is supplied with an Insertion Tip to help facilitate the insertion of the mesh pocket into the internal obturator muscle. See below for steps on how to use the Insertion Tip. Ensure that once the Sling is properly placed, the Insertion Tip is removed along with the forceps.

Contasure Needleless removal:

The removal technique of Contasure Needleless can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

- 1.- Vaginal incision: With the patient in modified dorsal lithotomy position, prepare the surgical field with a vaginal valve. Start by making a longitudinal incision in the anterior vaginal wall, preparing the submucosal space to present the best possible view of the urethral tract.
- 2.- Identify the mesh by direct view and palpation. Dissect the mesh from the surrounding structures. Avoid cutting the mesh whenever possible. Dissection of the mesh should continue by slightly pulling it to identify the mesh path. When mesh is completely dissected, remove the mesh from the vaginal side. Rinse the entire surgical site with antibiotic saline solution and close the vaginal incisions.

HOW TO USE THE INSERTION TIP**MATERIALS USED IN MANUFACTURING PROCESS**Patient-contacting materials:

Sling: Polypropylene monofilament mesh (colorless) and a blue polypropylene monofilament central traction thread (colorant: solvent blue 104, CAS no.: 116-75-6)
 Insertion Tip: Stainless steel* and UHMWPE (ultra high molecular weight polyethylene) multifilament

No phthalates, latex or other hazardous materials are used in the manufacture of Contasure Needleless mesh.

*Stainless steel components contain nickel

STERILIZATION

Contasure Needleless is sterilized by Ethylene Oxide. This product is a single use device and it must not be re-sterilized or re-used. Do not use if the packaging has been damaged.


STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and sunlight. Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT DISCLAIMER

Device intended to be used only once and for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON LABELING

 REF	CATALOGUE NUMBER		DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		KEEP AWAY FROM SUNLIGHT
 LOT	BATCH CODE		DISTRIBUTOR	 STERILE EO	STERILIZED USING ETHYLENE OXIDE
 SN	SERIAL NUMBER		CONSULT INSTRUCTIONS FOR USE		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE
 UDI	UNIQUE DEVICE IDENTIFIER		USE BY DATE	 31	DATE OF IMPLANTATION
 MD	MEDICAL DEVICE		MANUFACTURER		PATIENT NAME OR PATIENT ID
	DO NOT REUSE		DATE OF MANUFACTURE		NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER
	DO NOT RESTERILIZE		KEEP DRY		INFORMATION WEBSITE FOR PATIENTS

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