AGERS

INSTRUCCIONES DE USO
INSTRUCTIONS FOR USE

AGENESIS REMEEX SYSTEM

CE 0197

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax : +34-976-891872
E-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
E-mail: info@neomedic.com

Date of last revision: DECEMBER 2024

IMPORTANTE

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IMPORTANT

Please, read the Instructions for Use carefully.

Not following the indications of the user instructions may cause malfunctioning of the System or cause damage to the patient.

This document aims to explain the instructions for use of the Agenesis Remeex System. These instructions for use are not a substitute for training to correct vaginal agenesis. These instructions for use are not a document of reference to correct vaginal agenesis. Agenesis Remeex System should only be used by surgeons after undertaking specialized training, and with extensive experience in the treatment of pelvic floor diseases, more specifically in the use of this device, and following the instructions for use.

These instructions for use explain the common steps for implanting the system. Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may be also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.


INDICATIONS FOR USE

Agenesis Remeex System is a mechanical system to increase the length of the vaginal membrane in patients with vaginal agenesis.

DESCRIPTION OF THE DEVICE

The advantage of this system is that it permits to adjust the traction of the vaginal membrane step by step until the desired length is reached, as a way to avoid grafts or other complicated surgical techniques.

This device has been designed for a single patient use only. System components:

1.- Varitensor: Each day the manipulator top turn screw can be activated to gradually increase the length of the neovagina. The Varitensor is the part of the system that permits the Physician to tighten the Traction Thread by winding and continuous upward pressure of the Traction Ball on the vaginal vestibule which stretches the mucosa, leading to elongation of the cavity to a more normal depth. The Traction Thread is wound each day postoperatively until the neovagina is adequately formed (with each complete turn of the top turn screw, the Traction Thread length is wound 1mm on each side). The system has a retention mechanism that maintains the traction of the Traction Thread. The Traction Thread can only be loosened by turning the top turn screw counterclockwise.

2.- Base Plate: The Base Plate supports the Varitensor over the patient's abdomen. It also widens the Traction Thread angle to avoid pelvic structures.

3.- Soft support: It's a soft material support between the abdomen and the Base Plate.

4.- Traction thread: The Traction Thread connects the Traction Ball to the Varitensor.

5.- Traction Ball: It is a spherical piece that will form the neovagina by gradually lengthening the vaginal membrane.


6.- Diameter Maintenance Ball: To maintain the diameter of the neovagina that has been created by the Traction Ball.

7, 8, 9.- Agenesis needle: A ligature carrier type instrument specifically designed to help in passing the Traction Thread through the vaginal membrane. It is a canalized needle with a sheath and an internal guide.

10.- Vaginal Dilator set: The Vaginal Dilators are used by the patient postoperatively to maintain the neovagina space until regular sexual activity is established.

11.- Screwdriver: It is the element used to fix the Traction Thread to the Varitensor by tightening the front "Fixation Screw" on the Varitensor.

MECHANICAL OPERATION

ASSEMBLED AGENESIS NEEDLE

MECHANICAL FEATURES / SPECIFICATIONS

Varitensor function: Every turn of 360° of top turn screw winds the Traction Thread 1mm on each side.

Wind up capacity: The Varitensor can coil up more than 200mm of Traction Thread per each side.

Maximum traction stress: The Varitensor can support the force of 18 Newtons produced by the Traction Thread, with a minimum safety coefficient of 5.

CONTRAINDICATIONS

Do not use the device:

- In patients with pre-existing conditions that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- in patients who are undergoing anticoagulation treatments.
- in patients who suffer urinary tract infections or obstructions.
- in patients with renal insufficiency.
- in patients with immune diseases.
- If the patient shows signs of tissue necrosis.
- in patients who suffers cancer or have received radiation in the area of treatment.
- in patients with known sensitivity or allergy to device materials.

WARNINGS AND PRECAUTIONS

Pre-operative

- The patient should be counseled that there are alternatives to surgical treatment for vaginal agenesis, and the reasons for choosing surgical formation of neovagina should be explained. In appropriate cases, conservative and non-invasive treatment of self-dilation should be recommended before choosing surgical treatment.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications.
- Patient counseling should include a discussion that some complications associated with the procedure may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The physician must determine the risks that should be included in the patient informed consent.
- Users must be familiar with the "Vecchiotti" surgical technique and must be adequately trained in the placement of Agenesis Remeex System surgical procedure before employing the device.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- As the anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended location for temporary implant are planned and known for each individual patient. Employment of imaging methods before and after device placement may aid in proper placement and confirm absence of injury to non-target anatomical structures.
- The risks and benefits of performing the Agenesis Remeex System procedure in patients with the following conditions should be carefully considered by the physician:
 - Concomitant or previous pelvic floor surgery and, Patients with conditions that would compromise healing.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomy position.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or resterilize the single-use components of this device. All the Agenesis Remeex System components are single-use except the Vaginal Dilators. Vaginal Dilators are reusable by only one patient. Reject and do not use any unpackaged system.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Good surgical practices should be followed for management of contaminated or infected wounds.
- Cystoscopy should be performed to confirm bladder integrity after passage of the Traction Thread, as well as a digital rectal exam for detection of rectal perforation.
- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of the Traction Thread will minimize risks.
- The use of clamps could damage or weaken the strength of the Traction Thread. Damage to the Traction Thread could cause treatment failure.
- Do not turn the Manipulator without complete control. It is crucial to record the number of turns, the direction of the turns (both clockwise and counterclockwise) and the date that it has been performed.
- No part of the device system is a permanent implant.
- The discarded device components may be a potential biohazard. The Agenesis needle components have sharp tips, which could potentially protrude when not in use. Handle and dispose of the Agenesis needle, as well as the other discarded components in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.

Post-operative

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths, urinary incontinence, voiding dysfunction, problems defecating or any other abnormal symptoms to the surgeon as soon as possible.
- During the days after initial surgery, use proper techniques and conditions to avoid contamination and infection at the suprapubic site of Traction Thread withdrawal. If the gauze between the Soft Support and the abdomen needs to be replaced, cleaning must be performed by appropriate medical personnel.

- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Inform the patient that until regular penetrative sexual activity is established, or during periods of no penetrative sexual activity, a vaginal dilator should be worn at night to reinforce the neovagina and prevent contraction.
- It is recommended that patients avoid penetrative sexual intercourse until all surgical wounds have healed. As this may vary from patient to patient, the Physician shall determine when sexual relations may begin.

POTENTIAL ADVERSE REACTIONS

Risks associated to laparoscopic surgery for vaginal agenesis include:

- Acute inflammatory tissue reaction/swelling and transitory local irritation, allergic reaction or other adverse reactions to foreign bodies may occur.
- Infection, fever, inflammation/swelling or edema, discomfort and/or pain, adhesion formation, dysuria, bleeding, hematuria, fistulas, hematomas, hemorrhage, abscess, and/or seromas.
- As a foreign body, the device could stimulate a pre-existing infection.
- Urogynecologic problems related to micturition. These may include urethral obstruction, voiding dysfunction/de novo incontinence in the form of urgency, frequency or nocturia.
- Altered bladder and/or bowel function. Problems defecating.
- Mild or acute pain. Discomfort.
- Nerve damage, neuromuscular problems.
- Vaginal spotting, bleeding or abnormal discharge. These symptoms may be associated to wound site complications.
- Dyspareunia.
- Vaginal scarring.

As with all surgical interventions, the procedure may not be successful. Furthermore, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, removal of the device.

SURGICAL TECHNIQUE

Before using Agenesis Remeex System the doctor must read and understand the instructions for use.

Through laparoscopy technique, the Traction Thread, which has been previously tied to the traction ball, will be passed from the vaginal membrane, through the vesicorectal space into the abdominal cavity, through the extraperitoneal space and will ultimately traverse the anterior abdominal wall. Two suprapubic reference marks will be made to identify the exit marks of the Traction Thread. Once the Traction Thread is withdrawn suprapubically, it will pass through the holes of the Soft Support and the Base Plate and will be inserted into the Varitensor for their subsequent fixation.

Anesthetics and antibiotic therapy:

The surgical procedure to implant the Agenesis Remeex System can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

Surgical technique:

The device is made of a Traction Ball with a diameter of 25 mm (it has two holes through which the Traction Thread is inserted) and the Varitensor (to wind the Traction Thread). The Varitensor rests on the Base Plate (with holes where the Traction Thread works as a pulley). The Base Plate is placed suprapubically on abdominal wall of the patient.



The Traction Ball is placed between the urethra and the perineum, just at the level retro hymeneal recess. The PVDF suture (Traction Thread) attached.

In order to achieve extraperitoneal withdrawal of the Traction Thread, the use of a conventional laparoscope is required. A 10 mm intraumbilical trocar is for the scope insertion, while 2 trocars of 5 mm (placed at the anterior-superior iliac spine of each side of the pelvis) will permit the insertion of the instruments that will allow for the transversal cutting of the peritoneum. Once the pre-peritoneal space has been traversed, the bladder and the rectum are identified. Then a space is created at this level with the help of the assistant's finger which is pressing the vaginal membrane in order to identify the dissection zone.

After tunneling the dissection zone, the Agenesis Needle (ligature carrier) is used to insert the Traction Thread (which is secured to the Traction Ball), to the abdominal cavity from the vaginal membrane. The Traction Thread is passed to the abdominal cavity. Next, the 5mm shaft is removed until it disappears from the peritoneal hole and, without placing again the trocar; it is descended by the preperitoneal space to the dissected space. This procedure is not difficult if the hydro dissection is used and if the peritoneum is been detached and the trocar shaft has been moved further away from the underlying blood vessels. When the shaft appears in the central zone, where the peritoneum has been dissected, then a grasping forceps is inserted to grasp the Traction Thread that will be withdrawn through the suprapubic incisions, through which the laparoscope trocar was inserted. The Traction Thread on each side is inserted through the holes of the Soft Support and the Base Plate, acting as a pulley to establish connection to the Varitensor. The Soft Support should be placed below the Base Plate.

The extraperitoneal Varitensor has an adjustable top turn screw that through winding, creates continuous upward pressure of the Traction Ball on the vaginal vestibule which stretches the mucosa, leading to elongation of the cavity to a more normal depth. The Traction Thread, once in tension, may be wound upward daily in the postoperative period to achieve the neovagina.

The incisions are closed. It is recommended to apply Povidone-Iodine and medical gauze between the incisions and the Soft Support. Four days of catheterization will be required until the Traction Ball moves upwards 4cm inside the vagina. This is an indicator that the urethra won't be compressed.

Post-op traction regulation:

The following days after device implantation, the top turn screw of the Varitensor must be smoothly turned clockwise in order to wind the Traction Thread to create pressure on the vaginal membrane and begin to form the neovagina. The lengthening of the vaginal membrane will be gradual. The patient may remain in the hospital for the complete duration of the process. Alternatively and if the Physician(s) deems appropriate, the patient may be released from the hospital and continue the lengthening process under the Physician's indications. The tightening process will be determined by the medical team.

Every day the Traction Thread can be wound approximately 1 cm; to achieve this, the top turn screw of the Varitensor must undergo 3 complete turns.


Device Removal:

After 7-10 days the Traction Ball, will be approximately 7-10 cm inside the vulvar zone. Once confirmed, the Tracion Thread is cut extraperitoneally, thus by pulling downward on the bottom of the Traction Ball, it may be removed transvaginally. The vagina should be approximately the width of a wide finger, and have a depth of 7-10 cm. This procedure must be carried out in aseptic conditions to avoid infections.

Once the device has been removed, the patient must be instructed to wear a vaginal dilator of adequate size for several months at night. It is recommended that during the months following surgery, the dilator size is increased in thickness by 0.5 cm every month until reaching 3± cm.

Complete healing will be achieved in the months following surgery. Thus, it is recommended to lubricate the dilators with an antibiotic and healing ointment, as well as perform iodine washing once the dilators is withdrawn in the morning. The post-operative period should essentially be painless, without the need of painkillers. Furthermore, micturition should be normal once the catheter is removed from the urethra.

Penetrative sexual relations may begin when complete healing has been achieved and determined by the Physician.

Post-op Dilaton:

Postoperative dilation should be performed by the patient. Dilators are used to enhance the chance of success by maintaining the length of the neovagina. The Physician(s) may recommend gradually increasing the size (width) of the Dilators. The Dilators should be worn at night or for 8-10 hours daily during the first months following surgery, and until regular penetrative sexual activity is established. If continuous sexual activity is not maintained, a vaginal dilator should be worn at night to reinforce the neovagina and prevent contraction.

After daily use, the Dilator should be washed and sterilized with an antiseptic solution before re-use.

Medium and Long term Follow-Up:

Patient follow-up is advised, so that the Physician may analyze the evolution of the neovagina and to correct any possible abnormalities.

MATERIALS USED IN MANUFACTURING

Patient-contacting materials

Traction Ball and Diameter Maintenance Balls: Ultra-high molecular weight polyethylene (UHMWPE)

Vaginal dilators: Ultra-high molecular weight polyethylene (UHMWPE)

Traction Threads: Polyvinylidene fluoride (PVDF) monofilament

Agenesis needle: stainless steel*, polytetrafluoroethylene (PTFE)

No phthalates, latex or other hazardous materials are used in the manufacture of the Agenesis Remeex System patient-contacting components.

*Stainless steel components contain nickel

STERILIZATION

Agenesis Remeex System is sterilized by Ethylene Oxide.

All the Agenesis Remeex System components are single-use except the Vaginal Dilators. Vaginal Dilators are reusable by only one patient. The single-use components of the device must not be re-sterilized or re-used.

Do not use if the package has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light.

Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT

All the Agenesis Remeex System components, except the Vaginal Dilators, are intended to be used only once and for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER		DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE
	BATCH CODE		USE BY DATE
	SERIAL NUMBER		MANUFACTURER
	DO NOT REUSE		KEEP DRY
	DO NOT RESTERILIZE		KEEP AWAY FROM SUNLIGHT
	CONSULT INSTRUCTIONS FOR USE		STERILIZED USING ETHYLENE OXIDE

DIMA