SURGICAL MESH

- **NEOMESH SOFT**
- **NEOMESH STANDARD**

INSTRUCCIONES DE USO
INSTRUCTIONS FOR USE

REMESH SURGICAL MESH

CE 0197

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax : +34-976-891872
e-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
e-mail: info@neomedic.com

Date of last revision: November 2024

IMPORTANTE

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IMPORTANT

Please, read the Instructions for Use carefully.

Not following the indications of the user instructions may cause system malfunction or damage to the patient.

This document intends to explain the instructions for use of the REMESH SURGICAL MESH. Remesh Surgical Mesh should only be used by doctors instructed in the treatment of this type of pathology and more specifically in the use of this device, following the instructions for use.

These instructions for use explain the common steps to implant the system. The procedure could change when using the system with special techniques.

INTENDED USE

Remesh Surgical Mesh is intended for the treatment of fascial defects or damaged connective tissues where additional reinforcement or reconstruction is required to prevent or repair hernias.

DESCRIPTION OF THE DEVICE

Remesh Surgical Mesh is a knitted, non-absorbable polypropylene monofilament mesh.

This device has been designed for single use only.

WARNING AND PRECAUTIONS

Pre-operative:

- The device must only be used by doctors trained in the placement of the Remesh Surgical Mesh. It is mandatory to know the mechanism of action of the device.
- The device must only be implanted by doctors who are qualified to carry out the required surgical techniques.
- Patient counseling should include a discussion that there may be alternative treatments to non-absorbable mesh implant, that the mesh to be implanted is a permanent implant, and that some complications associated with the implanted mesh may require additional surgery that may or may not correct the complication.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- The physician must determine the risks that should be included in the patient informed consent.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to verify that it is not damaged and that the sterility is not compromised.
- This product is not reusable. Remesh Surgical Mesh cannot be re-sterilized. Reject and do not use any unpackaged system.
- The date of the last update of this document is included at the end of the Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant; therefore strict aseptic measures must be taken during the surgical procedure.
- Appropriate surgical procedures must be performed to implant Remesh Surgical Mesh as well as for the treatment of contaminated or infected wounds.
- Take special care during surgical intervention and the mesh fixation to avoid damage of internal structures, nerves, vessels, or mechanical disruption of the tissue and/or mesh material.
- Avoid direct bonding with internal organs, including viscera, to minimize the possibility of adhesions.
- Avoid excess tension of the mesh to minimize bowel obstruction or the possibility of damage to tissues or internal organs.

Post-operative:

- Advise the patient to contact the surgeon if any alterations or abnormal symptoms are detected.
- It is recommended that patients do not overly exert themselves or participate in strenuous physical activities or sports (for example biking, jogging, etc.) for a minimum of two months following surgery.

CONTRAINDICATIONS

Device is contraindicated:

- in patients with pre-existing conditions or in areas of the body that pose an unacceptable surgical risk. These conditions may include but are not limited to: patients with autoimmune diseases or renal insufficiency, patients undergoing anticoagulation treatments, patients suffering from infections, patients with known sensitivities or allergies to polypropylene, and in areas of the body affected by necrosis, cancer or radiation treatment.
- in the abdominal area in pregnant women or women that are considering future pregnancies.
- in babies or children, when its use may affect future growth.

ADVERSE REACTIONS

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur. These reactions may produce infection, inflammation, recurrence of the soft tissue defect, pain, fistulas, hematomas, abscess, erosion, extrusion and/or seromas.
- As a foreign body, Remesh Surgical Mesh could stimulate a pre-existing infection.
- Mesh contraction (shrinkage).
- Remesh Surgical Mesh may create adhesions when placed in direct contact with internal organs.
- Remesh Surgical Mesh may migrate if the mesh is not properly positioned or fixated.
- Foreign body sensitivity or palpability, discomfort.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to antibiotics, blood transfusions, intravenous therapy (IV), drainage, multiple reconstructive surgeries, and removal of the device.

SURGICAL TECHNIQUE

Before using Remesh Surgical Mesh the doctor must read and understand the instructions for use.

Anesthetics and antibiotic therapy:

The surgical procedure to implant the Remesh Surgical Mesh can be performed under local, regional or general anesthesia. It is advisable to follow the hospital antibiotic protocol for implants.

Surgical technique:

Rectangular/square mesh: mesh can be cut to a suitable size considering the anatomical location of the defect, the position of the patient, and his or her weight. Careful attention to mesh positioning and placement of fixation method will ensure optimal results. Fixation method must use only non-absorbable components. It is recommended that the fixation components are placed a minimum of 6mm from the edge of the mesh. Adequate spacing between fixation components will help prevent excessive tension or disruption between mesh material and connective tissue.

Preshape mesh: Careful attention to mesh shape and size selection, mesh positioning and placement of fixation method will ensure optimal results. The anatomical location of the defect, the position of the patient, and his or her weight must all be considered when selecting a mesh. Fixation method must use only non-absorbable components. It is recommended that the fixation components are placed a minimum of 6mm from the edge of the mesh. Adequate spacing between fixation components will help prevent excessive tension or disruption between mesh material and connective tissue.

Plug: Careful attention to plug size selection and positioning will ensure optimal results. Do not cut or trim the mesh plug as it is sized to fit and form to the irregular shape of the defect. Insertion of the mesh plug is accomplished by 1) grasping the center bulk portion of mesh with clean, sterile forceps, 2) centering the tapered end over the defect and, 3) gently pushing the plug into the hernia defect to the desired depth and position. Once inserted, the mesh plug should then be secured in position by fixating one of the knitted rows of the mesh plug to tissue. A flat polypropylene mesh onlay is provided with select models for placement directly over the inserted plug. When using a flat onlay mesh, it is recommended that the fixation components are placed a minimum of 6mm from the edge of the mesh. Fixation method must use only non-absorbable components. Adequate spacing between fixation components will help prevent excessive tension or disruption between mesh material and connective tissue.

MATERIALS USED IN MANUFACTURING

Polypropylene monofilament

Remesh Surgical Mesh is free of phthalates, latex and other hazardous materials.

STERILIZATION

Remesh Surgical Mesh is sterilized by Ethylene Oxide.
This product is a single use device and must not be re-sterilized or re-used.
Do not use if the packaging has been damaged.

STORAGE

Store the Remesh Surgical Mesh in a dry place at room temperature, away from heat and direct light. Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT DISCLAIMER

Device intended to be used once only for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER		CONSULT INSTRUCTIONS FOR USE		STERILIZED USING ETHYLENE OXIDE
	BATCH CODE		USE BY DATE		DATE OF IMPLANTATION
	SERIAL NUMBER		MANUFACTURER		PATIENT NAME OR PATIENT ID
	DO NOT REUSE		KEEP DRY		NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER
	DO NOT RESTERILIZE		KEEP AWAY FROM SUNLIGHT		INFORMATION WEBSITE FOR PATIENTS
	DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		DEVICE NAME		