DIMA

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IMPORTANT

Please, read the all relevant Instructions for Use carefully.
Not following user instructions or indications may cause malfunctioning of the device or cause damage to the patient.

This document aims to explain the instructions for use of DIMA Surgical Instruments for surgical placement of DIMA pelvic floor implant devices. These instructions for use are not a substitute for training to correct pelvic floor disorders. DIMA Surgical Instruments should only be used by doctors instructed in the treatment of pelvic floor disorders and more specifically in the use of compatible DIMA implant devices. User must consult the implant device's instructions for use before surgery.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective use of the device. As with any surgical intervention, special circumstances may be also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

INDICATIONS FOR USE

The device is designed to be used as surgical instrument during surgery, which aids in the insertion, placement and fixation of implant devices to correct disorders of the pelvic floor.

DESCRIPTION OF THE DEVICE AND DEVICE ELEMENTS

All components of this device have been designed for single use only. DIMA Surgical Instruments must only be used with compatible DIMA pelvic floor implant devices.

Commercial Ref.	Content & Description
APLIC03-D	Instrument for sacrospinous ligament fixation

The product reference listed above is intended to be used with the following compatible DIMA implant devices: Surelift Anchorsure, Surelift Contasure Prolapse System, Surelift Minimally Invasive Prolapse System.

CONTRAINDICATIONS

Do not use the product:

- in patients with known sensitivity or allergy to device materials.
- in patients contraindicated according to the instructions for use of the medical device together with which the surgical instrument is to be used.

WARNINGS AND PRECAUTIONS

Warnings and precautions regarding pelvic floor surgery are addressed in the implant device's instructions for use.

Pre-operative:

- The patient should be counseled about global mesh warnings before choosing surgical treatment for pelvic floor disorders.
- Users should be familiar with the surgical technique chosen, and should be adequately trained in the placement of the pelvic floor device implantation procedure before employing the device.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or reesterilize this device. No component of the DIMA Surgical Instrument is reusable. Reject and do not use any unpackaged system.
- The Anchoring Handle is supplied by the manufacturer without anchor/suture implant devices. The user must properly load the Anchoring Handle before use.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is used in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Cystoscopy should be performed to confirm bladder integrity, as well as a digital rectal exam for detection of rectal perforation.
- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of Anchoring Handle will minimize risks.
- Used devices may be a potential biohazard. DIMA Surgical Instruments could potentially protrude when not in use. Handle and dispose of DIMA Surgical Instruments, as well as the discarded components of the implant device in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Instruct the patient to report or any other abnormal symptoms (e.g.: heavy bleeding, intense pains) to the surgeon as soon as possible.

POTENTIAL ADVERSE REACTIONS

General risks related to the use of pelvic floor surgical instruments include:

- Acute inflammatory tissue reaction/swelling and transitory local irritation. Seroma.
- Nerve damage, neuromuscular problems; acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area. Dysuria. Other pains.
- Infection, stimulation of a pre-existing infection.
- Bleeding (hemorrhage, perforation of vessels, hematuria, edema, erythema, hematoma).

As with all surgical intervention, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence), scarring, physical deformity or emotional problems.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage.

SURGICAL USE

Before using DIMA Surgical Instruments, the doctor must read and understand the instructions for use included in the compatible DIMA implant device packaging.

Loading the Anchoring Handle:


The Anchoring Handle is composed of two parts: an inner rod and an outer tube, each with a polymer grip/trigger button. Inner: trigger button with two linear grooves; Outer: grip with round, central aperture.

The Anchoring Handle is designed specifically for a surgical technique called Sacrospinous Ligament Fixation (SSLF) and to be used with a SSLF system designed and manufactured by DIMA, composed of a standard surgical suture which is connected to a fixation anchor.

To load additional Anchors + Sutures:

1. Place the Anchor into the tip of the Anchoring Handle's inner rod, with each end of the suture passing through the lateral grooves in the tip.

2. With the Anchor loaded onto the tip of the Anchoring Handle's inner rod, insert the inner rod into the outer tube through the round opening in the grip.

3. Once the Anchoring Handle is assembled and loaded with an Anchor, be sure to insert the each end of the monofilament suture into the lateral grooves of the trigger button.

MATERIALS USED IN MANUFACTURING PROCESS

Anchoring Handle inner rod and outer tube: Stainless steel*
Anchoring Handle polymer grips: polypropylene, polyoxymethylene (POM)

**Stainless steel components contain nickel*

STERILIZATION

DIMA Surgical Instruments are sterilized with Ethylene Oxide.
This product is a single use device and must not be re-sterilized or re-used.
Do not use if the packaging has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light.
Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT

Device intended to be used once only for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER
	SERIAL NUMBER
	DO NOT REUSE
	DO NOT RESTERILIZE
	USE BY DATE
	STERILIZED USING ETHYLENE OXIDE

	DO NOT USE IF PACKAGE IS DAMAGED
	CONSULT INSTRUCTIONS FOR USE
	MANUFACTURER
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	SINGLE STERILE BARRIER SYSTEM

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