SURELIFT ANCHORSURE

**INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MODE D'EMPLOI
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
INSTRUCȚIUNI DE UTILIZARE
INSTRUÇÕES DE UTILIZAÇÃO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ**

NEOMEDIC INTERNATIONAL SL
C/ Maestrat 41-43, 08225 Terrassa (Barcelona), Spain
+34-93-7804505, info@neomedic.com, www.neomedic.com

DESARROLLO E INVESTIGACIÓN MÉDICA ARAGONESA S.L. (DIMA SL)
Polígono Industrial Mediavega, parcela 2.8, 50300 Calatayud (Zaragoza), Spain
+34-976-884065, dimasl@dimasl.com, www.dimasl.com

SURELIFT ANCHORSURE

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MODE D'EMPLOI
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
INSTRUCȚIUNI DE UTILIZARE
INSTRUÇÕES DE UTILIZAÇÃO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

IMPORTANT

Please read the Instructions for Use carefully.

Not following user instructions or indications may cause malfunctioning of the system or damage to the patient.

These instructions for use explain the common steps for implanting the Surelift Anchorsure. These instructions for use are not a substitute for training to correct pelvic organ prolapse (POP). The device should only be used by doctors instructed in the treatment of this condition, more specifically in the use of this device, and by following the instructions for use.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

The device summary of safety and clinical performance (SSCP) is made publically available through the European Database on Medical Devices (Eudamed), and provides healthcare workers and relevant patients critical information regarding the medical device in order to enhance the reader's understanding of the device's clinical safety and performance. The Eudamed public website can be accessed at the following address: <https://ec.europa.eu/tools/eudamed>. The SSCP for the device Surelift Anchorsure may be found by performing a search of its Basic UDI-DI: 843656835ANCHORSURENW.

Until the public Eudamed database module allowing access to the SSCP is up and running, the SSCP may be available to the public upon request, by contacting the manufacturer.

GLOBAL MESH WARNING

Surelift Anchorsure may be used together with a synthetic prolapse mesh, a device which is currently under heavy surveillance and monitoring. Health and governing authorities support the use of these devices when the mesh complies with certain specifications of synthetic meshes suitable for use as a pelvic floor implant (defined by the state of the art), when a multi-disciplinary team of medical experts considers surgical repair of pelvic organ prolapse (POP) with a mesh as an appropriate treatment option for the patient, and after considering all relevant risk factors. Surelift Anchorsure complies with the state of the art requirements for sacrospinous ligament fixation (SSLF). If Surelift Anchorsure is used together with a prolapse mesh implant, the mesh shall comply with the state of the art regarding physical and mechanical properties. Measures are in place globally in order to maximize the safety of prolapse mesh use, improve the availability of well-trained Physicians who can help patients decide if or when to use a mesh, and to more skillfully address the unintended outcomes associated with their use, with a focus on anatomic issues and ultimately demonstrating the benefit/risk profile of these devices at long-term.

INDICATIONS FOR USE

Surelift Anchorsure is indicated for attaching sutures to the sacrospinous ligament for vaginal prolapse repair.

DESCRIPTION OF DEVICE

Surelift Anchorsure is a sacrospinous-ligament (SSL) fixation system for pelvic organ prolapse treatment. It is composed of an anchor which is secured to the SSL using a specially designed applicator (Anchoring Handle), and an attached monofilament suture. The monofilament can be directly sutured to the prolapsed organ, or used to hold a synthetic mesh in place, creating further structural support, and successfully restoring the affected organs to their natural positions. Medical experts and health authorities recommend that synthetic prolapse meshes only be used for the treatment of moderate to severe, and/or complex cases of vaginal prolapse of the anterior or posterior compartment.

CONTRAINDICATIONS

Do not implant the device:

- in patients with pre-existing conditions that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- in patients who are undergoing anticoagulation treatments.
- in patients who suffer urinary tract infections or obstructions.
- in pregnant women or women that are considering future pregnancies.
- in patients with renal insufficiency.
- in patients with immune diseases.
- in patients who are still growing or less than 18 years of age.
- in areas with signs of tissue necrosis.
- in patients with known sensitivity or allergy to device materials.
- in patients undergoing concomitant bowel surgery.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to prolapse surgery and the treatment choice should be explained. Such alternatives include but are not limited to: Pelvic Floor Physical Therapy (PFPT; e.g.: Kegel exercises, electrostimulation), pessary/vaginal tampon, and or medication.
- Conservative, non-invasive treatments should be recommended before choosing surgical treatment. In cases of light or non-bothersome prolapse, conservative treatments should be exhausted before considering surgical repair. Except in cases of severe or extremely complex prolapse, special precautions should be taken for healthy patients under the age of 50; conservative treatments and/or native-tissue repair should be first-line treatment.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of prolapse surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implant may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The physician must determine the risks that should be included in the patient informed consent.

- Users should be familiar with the surgical technique(s) for POP treatment and should be adequately trained in the placement of Surelift Anchorsure implantation procedure before employing the device.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- As anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended location for product placement are planned and known for each individual patient. Employment of imaging methods before and after product placement may aid in proper product placement and confirm absence of injury to non-targeted anatomical structures.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomy position.
- The risks and benefits of performing the Surelift Anchorsure procedure in patients with the following conditions should be carefully considered by the team of Physicians:
 - Overweight patients;
 - Younger / pre-menopausal women;
 - Patients with conditions that would compromise healing, including but not limited to patients: who suffer cancers of the vagina, cervix, or uterus, who have received radiation in the area of treatment, with diabetes mellitus, who smoke or undergo steroid treatments;
 - Patients with $\leq$ grade 2 prolapse (POP-Q scale) considering a mesh;
 - Previous surgical history, especially vaginal surgeries for pelvic floor repair and including patients that already have a system implanted for prolapse or urinary incontinence;
 - Patients with concomitant SUI or urinary disorders.
- The patient should be counseled that future pregnancies may negate the effects of surgical procedure and the patient may suffer prolapse recurrence.
- Mesh devices are not recommended for women in the first two years after their last pregnancy.
- This product is presented sterile and must remain sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or resterilize this device. No component of the Surelift Anchorsure system is reusable. Reuse of the device could lead to potentially serious consequences for the patient that may lead to medical intervention, permanent injury or death. Reject and do not use any unpackaged system.
- Both the suture sleeve and the grip pin of the Anchoring Handle are components that must be discarded prior to the procedure.
- The date of the last revision of this document is included in this Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of product infection.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of the Anchoring Handle and other surgical instruments (if used) will minimize risks.
- Intra-operative bleeding may occur.
- Omitting the vaginal incision to insert the Anchoring Handle may result in excessive force and could cause injury to the patient.
- Cystoscopy should be performed to confirm bladder integrity, as well as a digital rectal exam for detection of rectal perforation.
- Mechanical damage may occur if affixing any part of the product with any staples, clips, or clamps. Clamps should only be used on the suture ends. The use of clamps could damage or weaken the strength of the suture. Cut and discard the area of the suture that is damaged. Damage to the suture could cause system failure.
- When connecting the suture to the surgical needle, only the suture end should be inserted into the eye. Insertion of the suture into the eye of the surgical needle could possibly damage the suture and cause system failure. Therefore, if the suture end has been damaged during insertion, the damaged area should be cut and discarded.
- Avoid excess tension of the Suture to minimize the possibility of Suture damage.
- No data is available on the use of tissue adhesives and is not recommended.
- Take special care to avoid tighten the Suture by pulling only one end. The Suture could disengage from the Anchor and it will no longer be possible to reconnect the Suture to the Anchor.
- The Anchoring Handle and the Surgical Needle are not permanent implants.
- The used components may be a potential biohazard. The Anchoring Handle and Surgical Needle have sharp tips, which could potentially protrude when not in use. Handle and dispose of the Anchoring Handle, the Surgical Needle, excess suture, and any other used components in accordance with any local and federal law regarding medical waste, with special attention to disposal requirement for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths, or any other abnormal symptoms to surgeon as soon as possible.
- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- It is recommended that patients do not overly exert themselves physically (e.g.: heavy lifting) or partake in any moderate to strenuous sporting activity (e.g.: biking, jogging, etc.) for at least a minimum of two months following surgery. It is also recommended to avoid sexual intercourse during the first two months after surgery.
- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as administered to the patient with the implant card for easy identification of the implant.
- Taking into account the condition of each individual patient, the surgeon will plan for post-operative follow-up. Follow-up is recommended at short/mid-term and long-term if possible.
- The Physician(s) and patient should report any serious incident in relation to Surelift Anchorsure to the manufacturer and the national health authority having jurisdiction in their locale.

POTENTIAL ADVERSE REACTIONS

Potential risks of sacrospinous ligament fixation surgery may include:

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur.
- Inflammation/swelling, edema, erythema, adhesion formation, fistulas, and/or seromas.
- Bleeding (hemorrhage, perforation of vessels, hematuria, hematoma)
- Perforation of bladder, bowel, vaginal wall, urethra, ureter, intestine and other pelvic structures.
- As a foreign body, Surelift Anchorsure could stimulate a pre-existing infection, or suffer rejection.
- Erosion/extrusion/exposure of the suture.
- Migration, erosion/extrusion/exposure, shrinkage, adhesion formation or hardening of the mesh caused by fibrosis (if used) may occur.

Associate mesh contraction (shrinkage) with vaginal shortening, vaginal tightening and vaginal pain.

- Urinary problems such as voiding symptoms, retention, de novo stress incontinence, or de novo incontinence in the form of urgency/overactive bladder (OAB), frequency or nocturia.
- Failure of the procedure resulting in recurrence (continued or worsening prolapse).
- Altered bowel function such as obstruction or problems defecating.
- Infection, wound infection, urinary tract infection, abscess.
- Neuromuscular problems, nerve damage; acute or chronic pain in the buttocks, lower back, thigh/leg and/or pelvic area. Dysuria. Other pains. Discomfort.
- Dyspareunia. Sexual partners may experience discomfort, irritation and/or pain during sexual activity.
- Vaginal spotting, bleeding, abnormal discharge or scarring.
- Foreign body sensitivity or palpability.

As with all surgical interventions, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

In cases in which the device and/or the mesh (if used) needs to be removed in part or whole, significant dissection may be required.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, and removal of the device and/or the mesh (if used).

SURGICAL TECHNIQUE

Before using the device, the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical technique to implant the device can be performed under regional or general anesthesia.

It is recommended to follow the surgical implants antibiotic protocol dictated by the hospital.

Surelift Anchorsure sacrospinous ligament fixation:

1.- Create a longitudinal incision in the anterior vaginal wall, extending from the midurethra to the vaginal apex (Figure A).

2.- Dissect anterior vaginal wall utilizing a full-thickness technique laterally to the level of tendinous arch on either side and proximally to the level of the vaginal apex. Perform anterior colporrhaphy if needed. Enter the paravesical spaces bilaterally. Identify ischial spines and clear them from adipose tissue digitally. Continue blunt/digital dissection into pararectal spaces bilaterally. Clear coccygeus muscle from adipose tissue digitally and identify course of sacrospinous ligament (SSL).

Note: Both the suture sleeve and the grip pin of the Anchoring Handle are components that must be discarded prior to the procedure.

3.- Insert the Anchoring Handle into the pararectal space on one side (Figure B) and position the tip over the Sacrospinous ligament (SSL) about 1.5 cm medial to the ischial spine (Figure C). With the tip of the Anchoring Handle firmly against the SSL, press the top of the Anchoring Handle as far as it allows, releasing the anchor and fixing it to the Sacrospinous ligament (Figure D). Withdraw the Anchoring Handle from the patient (Figure E) and repeat the technique on contralateral side if necessary (Figure F).

Figure A.

Figure B.

Figure C.

Figure D.



Figure E.

Figure F.

4. The sutures may then be employed by suspending the prolapsed vagina by direct suturing to the vaginal apex (vaginal vault) or by lifting the prolapsed organs from below with the aid of a synthetic mesh. In this manner, the suture can be tied to the posterior points of a synthetic mesh, holding it in place.

Surelift Anchorsure device removal:

In the event that removal of Surelift Anchorsure is required, careful considerations should be taken:

- 1.- If removal is deemed necessary because of chronic/persistent or bothersome pain, it is recommended that first, all sutures which have been fixed to pelvic organ structures (vaginal vault or uterus), are cut and removed.
- 2.- If the cause of removal is due to nerve or vessel damage stemming from the sacrospinous ligament, it is possible that the fixation anchors need to be explanted. Each individual case should be evaluated by the surgical team. The decision to explant fixation anchors is ultimately at the surgeon's discretion after a thorough assessment of the patient's state of health, given that additional risks may arise from removal of the anchors.

Loading the Anchoring Handle


The Anchoring Handle is composed of two parts: an inner rod and an outer tube, each with a polymer grip/trigger button. Inner: trigger button with two linear grooves; Outer: grip with round, central aperture.

Surelift Anchorsure system kit includes an assembled Anchoring Handle, pre-loaded with one Anchor + Suture. To load additional Anchor +

Suture:

1. Place the Anchor into the tip of the Anchoring Handle's inner rod, with each end of the suture passing through the lateral grooves in the tip.

2. With the Anchor loaded onto the tip of the Anchoring Handle's inner rod, insert the inner rod into the outer tube through the round opening in the grip.

3. Once the Anchoring Handle is assembled and loaded with an Anchor, be sure to insert each end of the monofilament suture into the lateral grooves of the trigger button.

MATERIALS USED IN MANUFACTURING

Patient contacting materials:

Anchor: Polyetheretherketone (PEEK)

Suture: Blue Polyvinylidene fluoride (PVDF) monofilament

Anchoring Handle inner rod and outer tube: Stainless steel*

Surgical Needle: Stainless steel*

No phthalates, latex or other hazardous materials are used in the manufacture of the Surelift Anchorsure Anchors and/or Sutures.

*Stainless steel components contain nickel

MAGNETIC RESONANCE SAFETY INFORMATION

MR Safe: The permanent implant components (Anchors and Sutures) of Surelift Anchorsure pose no safety hazards in the MR environment and may be placed anywhere in the MR environment. Patients with MR Safe devices have no scanning restrictions.

STERILIZATION

Surelift Anchorsure is sterilized by Ethylene Oxide. This product is a single use device and must not be re-sterilized or re-used. Do not use if the packaging has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and sunlight. Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT DISCLAIMER

Device intended to be used once only for a single patient. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

	CATALOGUE NUMBER		DISTRIBUTOR		STERILIZED USING ETHYLENE OXIDE
	SERIAL NUMBER		CONSULT INSTRUCTIONS FOR USE		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE
	UNIQUE DEVICE IDENTIFIER		USE BY DATE		DATE OF IMPLANTATION
	MEDICAL DEVICE		MANUFACTURER		PATIENT NAME OR PATIENT ID
	DO NOT REUSE		DATE OF MANUFACTURE		NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER
	DO NOT RESTERILIZE		KEEP DRY		INFORMATION WEBSITE FOR PATIENTS
	DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		KEEP AWAY FROM SUNLIGHT		MR SAFE

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