PASSERS

RETROPUBIC / TRANSOBTURATOR

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INSTRUCTIONS FOR USE

MODE D'EMPLOI

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DIMA SURGICAL INSTRUMENTS

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax : +34-976-891872
E-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
E-mail: info@neomedic.com

DATE OF LAST REVISION: JULY 2025

IMPORTANTE

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IMPORTANT

Please, read the all relevant instructions for use carefully.
Not following user instructions or indications may cause malfunctioning of the device or cause damage to the patient.

This document aims to explain the instructions for use of DIMA Surgical Instruments for surgical placement of DIMA pelvic floor implant devices. These instructions for use are not a substitute for training to correct pelvic floor disorders. DIMA Surgical Instruments should only be used by doctors instructed in the treatment of pelvic floor disorders and more specifically in the use of compatible DIMA implant devices. User must consult the implant device's instructions for use before surgery.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective use of the device. As with any surgical intervention, special circumstances may be also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

INDICATIONS FOR USE

The device is designed to be used as surgical instrument during surgery, which aids in the insertion, placement and fixation of implant devices to correct disorders of the pelvic floor.

DESCRIPTION OF THE DEVICE AND DEVICE ELEMENTS

All components of this device have been designed for single use only. The product may be composed of transobturator or retropubic surgical instruments. DIMA Surgical Instruments must only be used with compatible DIMA pelvic floor implant devices.

	Commercial Ref.	Content & Description
Retropubic Passer handles	MPAS11-D	For use with Ø3.5mm needle
	MPAS16-D	For use with Ø3.5mm needle
	MPAS17-D	For use with Ø5 mm needle
	MPAS09-D	For use with Ø3.5mm needle
	AR-05-D	For use with Ø5 mm needle
Retropubic Passer needles	PAS21-D	Length: 230mm Ø3.5mm, with t-slot tip
	PAS21-O-D	Length: 230mm Ø3.5mm, with oval slot tip
	PAS22-D	Length: 280mm Ø3.5mm, with t-slot tip
	AR-03-D	Length: 230mm Ø5mm, with circular slot tip
	PAS05-D	Length: 280mm Ø5mm, with circular loop slot tip
Transobturator Passer kits	APAS405-D	OUTSIDE-IN Helical (large) needles, Ø3.17mm
	APAS407-D	OUTSIDE-IN Helical (flare) needles, Ø3.17mm
	APAS410-D	OUTSIDE-IN Helical needles, Ø3.17mm
	APAS406-D	INSIDE-OUT Halo needles, Ø3.17mm
Retropubic Passer kits	APAS412-D	Length: 230mm Ø3.5mm, with t-slot tip
	TV35-21-D	Length: 230mm Ø3.5mm, with t-slot tip
	TV35-21-O-D	Length: 230mm Ø3.5mm, with oval slot tip
	TV5-PAS05-D	Length: 280mm Ø5mm, with circular slot tip

The product references listed are intended to be used with the following compatible DIMA implant devices: Contasure KIM, Remeex System for Female Urinary Incontinence, Remeex System for Male Urinary Incontinence and Surelift Contasure Prolapse System.

CONTRAINDICATIONS

Do not use the device:

- in patients with known sensitivity or allergy to device materials.
- in patients contraindicated according to the instructions for use of the medical device together with which the surgical instrument is to be used.

WARNINGS AND PRECAUTIONS

Warnings and precautions regarding pelvic floor surgery are addressed in the implant device's instructions for use.

Pre-operative

- The patient should be counseled about global mesh warnings before choosing surgical treatment for pelvic floor disorders.
- Users should be familiar with the surgical technique chosen, and should be adequately trained in the placement of the pelvic floor device implantation procedure before employing the device.
- This product is provided sterile and must be sterile before its use. Before opening the product the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or reesterilize this device. No component of the DIMA Surgical Instruments is reusable. Reject and do not use any unpackaged system.
- Before surgery, remove and discard the protective tips on the Passers and, if applicable, the winged Trocar.
- The detachable retropubic Passers are supplied unassembled by the manufacturer. Retropubic needles must be securely attached to the handle by fastening the turn screw before use.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is used in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Cystoscopy should be performed to confirm bladder integrity after the Passer needle passage, as well as a digital rectal exam for detection of rectal perforation.

- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of needles will minimize risks.
- Used devices may be a potential biohazard. DIMA Surgical Instruments have sharp tips, which could potentially protrude when not in use. Handle and dispose of DIMA Surgical Instruments, as well as the discarded components of the implant device in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.
- Omitting the vaginal and skin incisions to insert the Passers may result in excessive force and could cause injury to the patient.

Post-operative:

- Instruct the patient to report or any other abnormal symptoms (e.g.: heavy bleeding, intense pains) to the surgeon as soon as possible.

POTENTIAL ADVERSE REACTIONS

General risks related to the use of pelvic floor surgical instruments include:

- Acute inflammatory tissue reaction/swelling and transitory local irritation. Seromas.
- Nerve damage, neuromuscular problems; acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area. Dysuria. Other pains.
- Infection, stimulation of a pre-existing infection.
- Bleeding (hemorrhage, perforation of vessels, hematuria, edema, erythema, hematoma).

As with all surgical intervention, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence), scarring, physical deformity or emotional problems.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage.

SURGICAL USE

Before using DIMA Surgical Instruments, the doctor must read and understand the instructions for use included in the compatible implant device packaging.

The suture loops of the implant device must be passed through the slots in the proximal tips of the Passer needles.

Connecting the implant device sutures to Passers with closed circle or oval slot tip:

For the transobturator INSIDE-OUT technique, the winged trocar may be used as guidance for consistent, safe passage through the dissected path from the vaginal canal until reaching the obturator foramen.

MATERIALS USED IN MANUFACTURING PROCESS

Passer needle and Trocar: Stainless steel*

Passer handle and handle elements: Polypropylene and/or stainless steel

**Stainless steel components contain nickel*

STERILIZATION

DIMA Surgical Instruments are sterilized with Ethylene Oxide.

This product is a single use device and must not be re-sterilized or re-used.

Do not use if the packaging has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light.

Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT

Device intended to be used once only for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

CATALOGUE NUMBER

DISTRIBUTOR

STERILIZED USING ETHYLENE OXIDE

SERIAL NUMBER

CONSULT INSTRUCTIONS FOR USE

SINGLE STERILE BARRIER SYSTEM

UNIQUE DEVICE IDENTIFIER

USE BY DATE

INFORMATION WEBSITE FOR PATIENTS

MEDICAL DEVICE

DATE OF MANUFACTURE

DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE

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