UPLIFT

INSTRUCCIONES DE USO
INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
INSTRUÇÕES DE USO

SURELIFT UPLIFT

CE 0197

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax: +34-976-891872
E-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
E-mail: info@neomedic.com

Date of last revision: FEBRUARY 2026

IMPORTANTE

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IMPORTANT

Please, read the Instructions for Use carefully.

Not following user instructions indications may cause malfunctioning of the system or damage to the patient.

This document aims to explain the instructions for use of Surelift Uplift. These instructions are not a document of reference to correct vaginal prolapse.

Surelift Uplift should only be used by surgeons after undertaking specialized training, and with extensive experience in the treatment of pelvic floor disorders, more specifically in the use of this device, and following the instructions for use.

These instructions for use explain the common steps for implanting Surelift Uplift. Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

GLOBAL MESH WARNING

Prolapse meshes, specifically transvaginal mesh devices, are currently under heavy surveillance and monitoring. Though Surelift Uplift is considered a synthetic prolapse mesh device, *it is not* placed transvaginally. In fact, the use of a mesh is considered state of the art for surgical treatment of vaginal vault prolapse, of which Surelift Uplift is intended to treat. Nevertheless and generally speaking, health and governing authorities support the general use of mesh devices, when the following are true: the mesh complies with certain specifications of synthetic meshes suitable for use as a pelvic floor implant (defined by the state of the art), when a multi-disciplinary team of medical experts considers surgical repair of pelvic organ prolapse (POP) with a mesh as an appropriate treatment option for the patient, and after considering all relevant risk factors. Surelift Uplift complies with the state of the art requirements for a prolapse mesh implant with regards to physical and mechanical properties. Measures are in place globally in order to maximize the safety of prolapse mesh use, improve the availability of well-trained physicians who can help patients decide if or when to use a mesh, and to more skillfully address the unintended outcomes associated with their use, with a focus on anatomic issues and ultimately demonstrating the acceptable benefit/risk profile of these devices at long-term.

INDICATIONS FOR USE

Surelift Uplift is indicated for use as bridging material for sacrocolposuspension and/or sacrocolpopexy (laparotomy, laparoscopic or robotic approach) where surgical treatment for vaginal vault prolapse is warranted.

DESCRIPTION OF DEVICE

Surelift Uplift is a polypropylene monofilament mesh.

Surelift Uplift Mesh is designed to be used as a bridging material for sacrocolposuspension/ sacrocolpopexy where the prolapsed vaginal vault is supported with a suspensory bridge fixed to the anterior longitudinal ligament of the sacrum. This procedure may be performed via an open transverse abdominal route, or via robotic or laparoscopic route.

According to the product reference, five fixation anchors and an anchor applicator (Anchoring Handle) may be included. The function of the anchors is to fix the mesh to the sacral promontory.

CONTRAINDICATIONS

Do not implant the device:

- In patients with pre-existing conditions that pose an unacceptable surgical risk as determined by the physician, including comorbidities that may compromise tissue healing.
- In patients who are undergoing anticoagulation treatments.
- In patients who suffer urinary tract infections or obstructions.
- In the patient with renal insufficiency.
- In patients with immune deficiency.
- In patients who are still growing or less than 18 years of age.
- In patients who suffers cancers of the vagina, cervix, or uterus or in patients who have received radiation in the area of treatment.
- In pregnant women or women that are considering future pregnancies.
- In areas with signs of tissue necrosis.
- In patients with known sensitivity or allergy to system materials.
- In patients undergoing concomitant bowel surgery.
- In patients that already have a mesh system implanted for apical prolapse.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to prolapse surgery and the treatment choice should be explained. Such alternatives include but are not limited to: Pelvic Floor Physical Therapy (PFPT; e.g.: Kegel exercises, electrostimulation), pessary/vaginal tampon, and or medication.
- Conservative, non-invasive treatments should be recommended before choosing a mesh implant. In cases of light or non-bothersome prolapse, conservative treatments should be exhausted before considering mesh repair.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implant may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The physician must determine the risks that should be included in the patient informed consent.
- The device must only be used by doctors trained in the placement of the device. It is mandatory to understand the operational principles of the device.
- The device must only be implanted by doctors who are qualified, and experienced to carry out the required surgical technique.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- As anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended place for product placement are planned and known for each individual patient. Employment of imaging methods before and after product placement may aid in proper product placement and confirm absence of injury to non-targeted anatomical structures.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomic position.
- The risks and benefits of performing the Surelift Uplift procedure in patients with the following conditions should be carefully considered by the team of Physicians:
 - Obese patients;
 - Comorbidities which are risk factors for impaired tissue healing, including but not limited to: diabetes mellitus, smoking and steroid treatments;
 - Patients with $\leq$ grade 2 prolapse (POP-Q scale) considering a mesh;
 - Previous surgical history, especially vaginal surgeries for pelvic floor repair and patients that already have a system implanted for prolapse or urinary incontinence.
 - Patients with concomitant SUI or urinary disorders.
- This product is presented sterile and must be sterile before its use. Before opening the product the package must be inspected to check it is not damaged and sterility is not compromised.
- This product is not reusable. Surelift Uplift cannot be re-sterilized. Reject and do not use any unpackaged system.

- Remove the Anchoring Handle grip pin before the operation, and before loading the Anchoring Handle.
- The date of the last revision of this document is included in this Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of product infection.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of the Anchoring Handle and other surgical instruments (if used) will minimize risks.
- Avoid excess tension of the mesh to minimize urinary retention, bowel obstruction or the possibility of damage to tissues or internal organs.
- Intraoperative bleeding may occur.
- Do not affix any part of the product with any staples, clips, or clamps as mechanical damage to the product may occur. No data is available on the use of tissue adhesives and is not recommended.
- The used components may be a potential biohazard. Handle and dispose of the Anchoring Handle, excess anchors, and any other used components in accordance with any local and federal law regarding medical waste, with special attention to disposal requirement for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Patient could suffer post-operational pain in the abdomen and surrounding areas (e.g.: back, sacrum). Excessive or chronic pain should be reported to the Physician.
- Instruct the patient to report dysuria, hematuria, paresthesia, hematemesis or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- It is recommended that patients do not overly exert themselves physically (i.e.: heavy lifting) or partake in any moderate to strenuous sporting activity (e.g.: biking, jogging, etc.) for at least a minimum of two months following surgery. It is also recommended to avoid sexual intercourse during the first two months after surgery.
- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as administered to the patient for easy identification of the implant.

POTENTIAL ADVERSE REACTIONS**Potential risks of prolapse mesh surgery include:**

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur. These reactions may produce fever, infection, inflammation, recurrence of the pathology, discomfort, temporary or chronic pain, dysuria, hematuria, fistulas, hematomas, hemorrhage, erythema, abscesses, and/or seromas. Hardening of the mesh caused by fibrosis.
- As a foreign body, Surelift Uplift could stimulate a pre-existing infection.
- Extrusion/exposure of the mesh or sutures, erosion of surrounding tissues (e.g.: vaginal).
- Migration or shrinkage of the mesh may occur. Associate mesh contraction (shrinkage) with vaginal shortening, vaginal tightening and vaginal pain.
- De novo or recurrent stress urinary incontinence, urge incontinence and/or urinary frequency may appear after the intervention.
- Failure of the procedure resulting in continued or worsening of the pathology. Recurrent or de novo prolapse of another compartment could occur.
- Altered bladder and/or bowel function. Bowel obstruction. Signs and symptoms of bowel obstruction include but are not limited to constipation, abdominal pain, vomiting and bloating.
- Neuromuscular problems, nerve damage.
- Dyspareunia.
- Redness, itching or dryness of the vagina. Vaginal discharge, bleeding.
- Foreign body sensitivity, discomfort.
- Vaginal scarring.

As with all surgical interventions, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence, infection, and herniation).
- Scarring.
- Physical deformity.
- Emotional problems.
- Post-op nausea and/ or vomiting.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, and removal of the mesh and/or anchors (if used). In cases in which the mesh and/or anchors (if used) needs to be removed in part or whole, significant dissection may be required.

SURGICAL TECHNIQUE

Before using Surelift Uplift the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical technique to implant Surelift Uplift can be performed under regional or general anaesthesia. It is recommended to follow the antibiotic protocol dictated by the hospital.

Surelift Uplift placement:

- 1.- If using laparoscopic approach, ensure compatibility by inserting the device into the trocar prior to the prolapse surgery. The Surelift Uplift Anchoring Handle should enter and exit the trocar easily without excessive force. The use of excessive force could damage the instrument.
- 2.- If the uterus is preserved, temporarily fix the uterus to the abdominal wall through percutaneous suture with a straight needle. Pass the suprapubic suture through the abdominal wall, passing through the uterine fundus and out again through the abdominal wall. For security, tie the suture onto a piece of gauze.
- 3.- Identify the sacral promontory and create an incision in the peritoneum, avoiding the nearby vessels.
- 4.- Continue the incision of the peritoneum lateral to the rectum to expose the pelvic muscles at the bottom of the pelvis. The peri-rectal fat must be preserved to avoid damage. Repeat dissection on the opposite side.
- 5.- Fix the Y-shaped end of the Uplift mesh bilaterally to the levator ani muscles. It is recommended to use a method for the correct identification of the vaginal dome (balloon probe).
- 6.- Pull the central part of the mesh up towards the promontory, applying the adequate tension and fix the mesh with anchors or sutures, depending on the surgeon's preference, to the vertebral anterior ligament.
- 7.- Remove the suture which temporarily holds the uterus to the abdominal wall and create an opening in the broad ligament of the uterus (right side), avoiding the uterine artery and Fallopian tube.
- 8.- Identify the top of the vaginal dome to guide the dissection of the bladder from the vagina. Dissect until reaching the balloon probe (in the event that this technique is being used).

9.- Pass the loose end of the mesh through the opening of the broad ligament, making sure that it is not wrapped up in itself, and that it remains this way during the placement. Fix it to the area of former dissection, the apex of the mesh being fixed or sutured on both sides.

10.- With forceps, pull the central part of the mesh towards the promontory once more, applying the adequate tension to correct the prolapse, and fix the mesh using two anchors or sutures to the anterior vertebral ligament.

11.- Cut the excess of mesh and close the peritoneum incision, avoiding mesh exposure.

12.- After use, the Anchoring Handle may be a potential biohazard. This device has a sharp tip, which could potentially protrude when not in use. Handle and dispose of the Anchoring Handle, as well as any unused or removed anchors in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices to prevent unnecessary injury or harm.




Loading the Anchoring Handle

The Anchoring Handle is composed of two parts: an inner rod and an outer tube, each with a polymer grip/trigger button.

Inner: colored trigger button; Outer: grip with round, central aperture.

Some Surelift Uplift product references include a preassembled Anchoring Handle (preassembled=the inner rod is inserted into the outer tube).

To load Uplift Anchors:

1.- The Anchoring Handle will be preassembled. Before use, remove the grip pin secured between the two Anchoring Handle grips. 	2.- With the assembled* Anchoring Handle, insert the Anchoring Handle outer tube into the Anchor carrier, ensuring that the Anchoring Handle outer tube is pushed down completely. *IMPORTANT: If assembled correctly, the sharp tip of the inner rod should not protrude from the outer tube (as seen in the third image marked by an "x"). 	3.- With the Anchoring Handle outer tube still inside the Anchor carrier and holding the Anchoring Handle by the larger grip, push slightly down on the top, smaller trigger button until halfway engaged inside the larger grip. IMPORTANT: The smaller trigger button should not be pushed all the way down (as seen in the third image marked by an "x"). 	4.- Pull the Anchoring Handle out from the carrier. The Anchor will stay fixed to the Anchoring Handle (the Anchor should not protrude from the Anchoring Handle outer tube). 
--	---	--	---

MATERIALS USED IN MANUFACTURING

Patient-contacting materials

Prolapse mesh: Polypropylene monofilament

Anchors: Polyetheretherketone (PEEK)

Anchoring Handle (tip): Stainless Steel*, ultra high molecular weight polyethylene (UHMWPE)

No phthalates, latex or other hazardous materials are used in the manufacture of Surelift Uplift mesh or anchors.

*Stainless steel components contain nickel

STERILIZATION

Surelift Uplift is sterilized by Ethylene Oxide.

This product is a single use device and it must not be re-sterilized or re-used. Do not use if the package has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light. Do not use after the expiration date detailed on the package.

SINGLE-USE PRODUCT DISCLAIMER

Device intended to be used only once and for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

CATALOGUE NUMBER

SERIAL NUMBER

UNIQUE DEVICE IDENTIFIER

MEDICAL DEVICE

DO NOT REUSE

DO NOT RESTERILIZE

DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE

DISTRIBUTOR

CONSULT INSTRUCTIONS FOR USE

USE BY DATE

MANUFACTURER

KEEP DRY

DATE OF MANUFACTURE

KEEP AWAY FROM SUNLIGHT

STERILIZED USING ETHYLENE OXIDE

SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE

DATE OF IMPLANTATION

PATIENT NAME OR PATIENT ID

NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER

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