
SURELIFT

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INSTRUCTIONS FOR USE

MODE D'EMPLOI

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SURELIFT CONTASURE PROLAPSE SYSTEM

CE 0197

MANUFACTURING

Desarrollo e Investigación Médica Aragonesa S.L.
Polígono Industrial Mediavega, Parcela 2.8
50300 Calatayud (Zaragoza)
Spain
Tel: +34-976-884065
Fax : +34-976-891872
e-mail: dimasl@dimasl.com

DISTRIBUTION

Productos Especiales de Neomedic International S.L.
C/ Maestrat, nº 41-43 1º
08225 Terrassa (Barcelona)
Spain
Tel: +34-93-7804505
Fax: +34-93-7806105
e-mail: info@neomedic.com

Date of last revision: JUNE 2025

IMPORTANTE

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ENGLISH

IMPORTANT

Please read the Instructions for Use carefully.

Not following the indications of the user instructions may cause malfunctioning of the system or damage to the patient.

This document aims to explain the instructions for use of the Surelift Contasure Prolapse System. These instructions are not a document of reference to correct the vaginal prolapse. Surelift Contasure Prolapse System should only be used by surgeons after undertaking specialized training, and with extensive experience in the treatment of pelvic floor disorders, more specifically in the use of this device, and following the instructions for use.

These instructions for use explain the common steps for implanting the system. Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

GLOBAL MESH WARNING

Surelift Contasure Prolapse System is considered a synthetic prolapse mesh device, of which are currently under heavy surveillance and monitoring. Health and governing authorities support the use of these devices when the following are true: the mesh complies with certain specifications of synthetic meshes suitable for their use as a pelvic floor implant (defined by the state of the art), when a multi-disciplinary team of medical experts considers surgical repair of pelvic organ prolapse (POP) with a mesh as an appropriate treatment option for the patient, and after considering all relevant risk factors. Surelift Contasure Prolapse System complies with the state of the art requirements for a prolapse mesh implant with regards to physical and mechanical properties. Measures are in place globally in order to maximize the safety of prolapse mesh use, improve the availability of well-trained Physicians who can help patients decide if or when to use a mesh, and to more skillfully address the unintended outcomes associated with their use, with a focus on anatomic issues and ultimately demonstrating the benefit/risk profile of these devices at long-term.

INDICATIONS FOR USE

Surelift Contasure Prolapse System is a mesh system for the treatment of moderate to severe, and/or complex cases of vaginal prolapse of the anterior or posterior compartment.

DESCRIPTION OF THE DEVICE

Surelift Contasure Prolapse System is a polypropylene monofilament mesh. It remains implanted reinforcing the pelvic floor by inducing a fibrotic tissue response through and around the mesh. The posterior anchoring points will be used to attach the mesh to the sacrospinous ligament (SSL). Anchors will be inserted into the SSL, functioning to secure the monofilament thread suture that will keep the posterior part of the mesh in place. The arms of the mesh will be passed through the obturator foramen.

CONTRAINDICATIONS

Do not implant the device:

- In patients with pre-existing conditions that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- In patients who are undergoing anticoagulation treatments.
- In patients who suffer urinary tract infections or obstructions.
- In the patient with renal insufficiency.
- In patients with immune deficiency.
- In patients who are still growing or less than 18 years of age.
- In patients who suffers cancers of the vagina, cervix, or uterus or in patients who have received radiation in the area of treatment.
- In pregnant women or women that are considering future pregnancies.
- In sub-urethral or peri-urethral areas with signs of tissue necrosis.
- In patients with known sensitivity or allergy to system materials.
- In patients undergoing concomitant bowel surgery.
- In direct contact with bowel or visceral organs, including the urinary bladder.
- In patients that already have a system implanted for prolapse or urinary incontinence.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to prolapse-mesh surgeries, and the reason for choosing a prolapse mesh should be explained. Such alternatives include but are not limited to: Pelvic Floor Exercises (Kegel exercises), Pessary/Vaginal Tampon, and/or Behavioral Modification (electrostimulation). The patient may also be offered non-mesh surgical options (e.g.: native-tissue repair).
- Conservative, non-invasive treatments should be recommended before choosing a mesh implant. In cases of light or non-bothersome prolapse, conservative treatments should be exhausted before considering mesh repair. Except in cases of severe or extremely complex prolapse, otherwise healthy patients under the age of 50 should only be offered conservative treatments and/or native-tissue repair.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implant may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The physician must determine the risks that should be included in the patient informed consent.
- The device must only be used by doctors trained in the placement of the device. It is mandatory to understand the operational principles of the device.
- The device must only be implanted by doctors who are qualified and experienced to carry out the required surgical technique.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- As anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended place for product placement are planned and known for each individual patient. Employment of imaging methods before and after product placement may aid in proper product placement and confirm absence of injury to non-targeted anatomical structures.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomic position.
- The risks and benefits of performing the Surelift Contasure Prolapse System procedure in patients with the following conditions should be carefully considered by the team of Physicians:
 - Obese patients;
 - Younger / pre-menopausal women or women planning future pregnancies;
 - Comorbidities which are risk factors for impaired tissue healing, including but not limited to: diabetes mellitus, smoking and steroid treatments;
 - Patients with $\leq$ grade 2 prolapse (POP-Q scale);
 - Previous surgical history, specially vaginal surgeries for pelvic floor repair;
 - Patients with concomitant SUI or urinary disorders.
- The patient should be counseled that future pregnancies may negate the effects of surgical procedure.
- This product is presented sterile and must be sterile before its use. Before opening the product the package must be inspected to check it is not damaged and sterility is not compromised.
- This product is not reusable. The system cannot be re-sterilized. Reject and do not use any unpackaged system.
- The date of the last revision of this document is included in this Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of product infection.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of the mesh will minimize risks.
- Cystoscopy should be performed to confirm bladder integrity after the needle-passage, as well as a digital rectal exam for detection of rectal perforation.
- Avoid excess tension of the mesh to minimize urinary retention, bowel obstruction or the possibility of damage to tissues or internal organs.
- Intraoperative bleeding may occur.
- Do not affix any part of the product with any staples, clips, or clamps as mechanical damage to the product may occur. No data is available on the use of Surelift Contasure Prolapse System with tissue adhesives and is not recommended.
- Improperly attaching the Suture to the surgical needle could cause damage to the Suture. If securing the Suture to the surgical needle through its open-ended slot, this should only be done with the Suture end.
- Cut and withdraw the mesh traction-threads before closing the incisions.
- The used components may be a potential biohazard. The Passers have sharp tips, which could potentially protrude when not in use. Handle and dispose of the Anchoring Handle, Passers, discarded mesh Traction Threads, and any other used components in accordance with any local and federal law regarding medical waste, with special attention to disposal requirement for sharp devices, to prevent unnecessary injury or harm.

SURELIFT CONTASURE PROLAPSE SYSTEM

INSTRUCTIONS FOR USE

ENGLISH

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative retropubic bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- It is recommended that patients do not overly exert themselves physically (e.g.: heavy lifting) or partake in any moderate to high intensity sporting activity (e.g.: biking, jogging, etc.) for at least a minimum of two months. It is also recommended that patients avoid sexual intercourse during the first two months after surgery.
- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as administered to the patient for easy identification of the implant.

POTENTIAL ADVERSE REACTIONS

Known risks of surgical procedures for prolapse mesh implanting include:

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur. These reactions may produce infection, inflammation, recurrence of the pathology, pain and/or discomfort, adhesion formation, dysuria, hematuria, fistulas, hematomas, hemorrhage, erythema, abscesses, and/or seromas.
- As a foreign body, the device could stimulate a pre-existing infection.
- Erosion, migration, extrusion or shrinkage of the mesh may occur.
- Calcification of device.
- De novo urge incontinence and/or urinary frequency may appear after the intervention.
- Failure of the procedure resulting in continued or worsening prolapse.
- Altered bladder and/or bowel function.
- Neuromuscular problems, nerve damage.
- Atrophic vaginitis which could be characterized by redness, itching or dryness of the vagina.
- Dyspareunia.
- Vaginal discharge, bleeding.
- Foreign body sensitivity or palpability, discomfort.

As with all surgical interventions, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, and removal of the device. In cases in which the polypropylene mesh needs to be removed in part or whole, significant dissection may be required.

SURGICAL TECHNIQUE

Before using Surelift Contasure Prolapse System, the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical procedure to implant the system can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

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Surelift Contasure Prolapse System placement: (Anterior Prolapse)

- 1.- Perform a longitudinal incision in the anterior vaginal wall extending from mid-urethra to the vaginal apex.
- 2.- Dissect the anterior vaginal wall utilizing a full-thickness technique laterally to the level of arcus tendineus on either side and proximally to the level of vaginal apex. Perform anterior colporrhaphy if needed. Enter paravesical spaces bilaterally. Identify ischial spines and clear them from adipose tissue digitally. Continue blunt/digital dissection into pararectal spaces bilaterally. Clear coccygeus muscle from adipose tissue digitally and identify course of sacrospinous ligament.
- 3.- Insert the Anchoring Handle into the pararectal space on one side and position the tip over the sacrospinous ligament about 1.5 cm medial to the ischial spine. With the tip of the Anchoring Handle firmly against the ligament, press the top of the Anchoring Handle as far as it allows, releasing the anchor and fixing it to the sacrospinous ligament. Withdraw the handle from the patient and secure the proximal traction sutures. Repeat the technique on contra lateral side.
- 4.- Put the suture through the one of the posterior anchoring points of the mesh (the shorter arms at the base of the mesh which do not have traction thread loops attached) and tie a "pulley" knot. Pull one end of the suture and bring the same arm of the mesh into direct contact with the desired ligament fixation point. Repeat the same steps on contra lateral side and adjust width of the mesh to patient's dimensions leaving it somewhat loose to accommodate for possible mesh shrinkage but avoiding excessive sagging and mesh folding. Tie the knot of the anchor threads and cut the excess suture ends.
- 5.- Perform four small vulvar incisions (2 on each side) lateral to labia majora. The first incisions should be approximately at the level of the clitoris and the medial-superior corner of the obturator foramen. The second incisions should be about 2 cm under and 2 cm lateral to the first incisions, inferior to the obturator foramen.
- 6.- Insert the thread Passer Needle through the vulvar incision and penetrate into the obturator membrane. Advance the tip of the Passer Needle through the vaginal canal. Insert and secure the suture in the apex of the tip of the Passer Needle. Pull the anterior arm of the mesh through the obturator foramen and anterior vulvar incision. Repeat the same steps on contra lateral side and adjust width of the mesh to patient's dimensions leaving it somewhat loose to accommodate for possible mesh shrinkage but avoiding excessive sagging and mesh folding. Adjustment can be facilitated by cutting any excess mesh in order to ensure the proper position of the mesh without tension under the bladder neck, leaving adequate space between the mesh and bladder neck to accommodate for possible mesh shrinkage.
- 7.- Insert the Passer Needle through the second vulvar incisions and penetrate until the obturator membrane is reached. Laterally advance the tip of the Passer Needle to the obturator internus and penetrate through arcus tendineus as close to ischial spine as possible. Advance the tip of the Passer Needle into the vaginal canal. Insert and secure the suture in the apex of the Passer Needle and pull the central arm of the mesh through obturator foramen and the second vulvar incision. Repeat the same steps on contra lateral side. Adjust width of the mesh to patient's dimensions leaving it somewhat loose to accommodate for possible mesh shrinkage but avoiding excessive sagging and mesh folding.
- 8.- Perform appropriate vaginal vault suspension as per surgeon's preference and patient's indications. Cut the mesh arms at skin level and close all incisions, including the anterior wall vaginal incision.

Surelift Contasure Prolapse System placement: (Posterior Prolapse)

- 1.- Perform a transverse incision in the posterior fourchette of the vagina or inverted triangle incision on the perineum as required by patients conditions.
- 2.- Dissect posterior vaginal wall sharply and bluntly, underlying the rectum, utilizing full-thickness technique all the way to vaginal apex proximally and to the anal elevator muscle fascia laterally, without making a midline incision. Continue blunt/digital dissection into pararectal spaces bilaterally. Clear coccygeus muscle from adipose tissue digitally and identify course of sacrospinous ligament.
- 3.- Insert the Anchoring Handle into the pararectal space on one side and position the tip over the Sacrospinous ligament (SSL) about 1.5 cm medial to the ischial spine. With the tip of the Anchoring Handle firmly against the SSL, press the top of the Anchoring Handle as far as it allows, releasing the anchor and fixing it to the Sacrospinous ligament. Withdraw the Anchoring Handle from the patient and secure the proximal sutures threads. Repeat the technique on contralateral side.
- 4.- Prepare the mesh, cutting the central arms and the anterior arms as per surgeon preference (they are not necessary for posterior placement). Put the suture through the one of the posterior anchoring points of the mesh (the shorter arms at the base of the mesh which do not have traction-threads) and tie a "pulley" knot. Pull the suture and bring the same arm of the mesh into direct contact with the desired ligament fixation point. Repeat the same steps on contra lateral side and adjust width of the mesh to patient's dimensions leaving it somewhat loose to accommodate for possible mesh shrinkage but avoiding excessive sagging and mesh folding. Tie the knot of the anchor threads to fix the mesh at the desired level, and cut the excess suture ends. Perform colpoperineorrhaphy if needed.
- 5.- Adjust the length of the mesh to patient's dimensions leaving it somewhat loose to accommodate for possible mesh shrinkage but avoiding excessive sagging and mesh folding. Suture anterior portion of the mesh to proximal portion of the perineal body with delayed absorbable sutures but do not extend mesh to posterior fourchette to prevent dyspareunia. Approximate posterior vaginal wall mucosa with perineal skin.
- 6.- Close all incisions, ensuring that only the mesh remains implanted.


SURELIFT CONTASURE PROLAPSE SYSTEM

INSTRUCTIONS FOR USE


ENGLISH

Loading the Anchoring Handle

The Anchoring Handle is composed of two parts: an inner rod and an outer tube, each with a polymer grip/trigger button. Inner: trigger button with two linear grooves; Outer: grip with round, central aperture. Surelift Contasure Prolapse System kits that include an assembled Anchoring Handle, are pre-loaded with one Anchor + Suture. To load additional Anchor + Suture:

1. Place the Anchor into the tip of the Anchoring Handle's inner rod, with each end of the suture passing through the lateral grooves in the tip.

2. With the Anchor loaded onto the tip of the Anchoring Handle's inner rod, insert the inner rod into the outer tube through the round opening in the grip.

3. Once the Anchoring Handle is assembled and loaded with an Anchor, be sure to insert the each end of the monofilament suture into the lateral grooves of the trigger button.

MATERIALS USED IN MANUFACTURING

Patient-contacting materials

Mesh: Polypropylene monofilament
Anchors: Polyetheretherketone (PEEK), polyvinylidene fluoride (PVDF) suture thread
Anchoring Handle*: Stainless Steel
Passers*: Stainless Steel
Surgical needle*: Stainless steel

No phthalates, latex or other hazardous materials are used in the manufacture of Surelift Contasure Prolapse System mesh, anchors or sutures.

*Stainless steel components contain nickel

STERILIZATION

Surelift Contasure Prolapse System is sterilized by Ethylene Oxide. This product is a single use device and it must not be re-sterilized or re-used. Do not use if the package has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and direct light. Do not use after the expiration date detailed on the package.

SINGLE-USE PRODUCT DISCLAIMER

Device intended to be used only once for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON PRODUCT LABELS

CATALOGUE NUMBER

SERIAL NUMBER

UNIQUE DEVICE IDENTIFIER

MEDICAL DEVICE

DO NOT REUSE

DO NOT RESTERILIZE

DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE

DISTRIBUTOR

CONSULT INSTRUCTIONS FOR USE

USE BY DATE

MANUFACTURER

KEEP DRY

DATE OF MANUFACTURE

KEEP AWAY FROM SUNLIGHT

STERILIZED USING ETHYLENE OXIDE

SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE

DATE OF IMPLANTATION

PATIENT NAME OR PATIENT ID

NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER

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