FEMALE REMEEX SYSTEM

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
GEBRAUCHSANWEISUNG
MODE D'EMPLOI
ISTRUZIONI PER L'USO
INSTRUÇÕES DE USO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

NEOMEDIC INTERNATIONAL SL
C/ Maestrat 41-43, 08225 Terrassa (Barcelona), Spain
+34-93-7804505, info@neomedic.com, www.neomedic.com

DIMA

DESARROLLO E INVESTIGACIÓN MÉDICA ARAGONESA S.L. (DIMA SL)
Polígono Industrial Mediavega, parcela 2.8, 50300 Calatayud (Zaragoza), Spain
+34-976-884065, dimasl@dimasl.com, www.dimasl.com

FEMALE REMEEX SYSTEM

INSTRUCCIONES DE USO
INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
MODE D'EMPLOI
ISTRUZIONI PER L'USO
INSTRUÇÕES DE USO
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

IMPORTANT

Please, read the Instructions for Use carefully.

Not following user instructions or indications may cause malfunctioning of the system or cause damage to the patient.

These instructions for use explain the common steps for implanting Remeex System for Female Urinary Incontinence. These instructions for use are not a substitute for training to correct SUI (Stress Urinary Incontinence). The device should only be used by doctors instructed in the treatment of this condition, more specifically in the use of this device, and by following the instructions for use.

Considering that each individual patient may present special anatomical conditions, it is imperative that the surgeon understands and plans for the safe and effective placement of the device. As with any surgical intervention, special circumstances may also be presented during the surgery; in such cases the medical specialist will be responsible for deviating from the instructions as necessary to ensure patient safety and health.

The device summary of safety and clinical performance (SSCP) is made publically available through the European Database on Medical Devices (Eudamed), and provides healthcare workers and relevant patients critical information regarding the medical device in order to enhance the reader's understanding of the device's clinical safety and performance. The Eudamed public website can be accessed at the following address: <https://ec.europa.eu/tools/eudamed>. The SSCP for the device Remeex System for Female Urinary Incontinence may be found by performing a search of its Basic UDI-DI: 843656835REMEEXFEMALEQN.

Until the public Eudamed database module allowing access to the SSCP is up and running, the SSCP may be available to the public upon request, by contacting the manufacturer.


INDICATIONS FOR USE

Remeex System for Female Urinary Incontinence is intended to be used as a suburethral sling for recurrent or moderate to severe female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

DESCRIPTION OF THE DEVICE AND DEVICE ELEMENTS

All components of this device have been designed for single use only. System components:

- 1.- Varitensor: The permanently implanted part of the system which permits the physician to adjust the Sling tension.
- 2.- Manipulator: The part of the system which permits the Varitensor to be adjusted from outside the body to increase or decrease the Sling tension. Before the patient is discharged from the hospital, the Manipulator is disconnected from the Varitensor and removed from the patient using the Disconnecter (the Manipulator comes connected to the Varitensor from the manufacturer).
- 3.- Disconnecter: The part of the system that is used to disconnect or connect the Manipulator to the Varitensor. The system includes 2 Disconnecters.
- 4.- Sling: A 15 x 32 mm polypropylene monofilament mesh connected to two monofilament Traction Threads. The Sling will be connected to the Varitensor by the Traction Threads.
- 5.-Traction Threads: The Sling will be connected to the Varitensor by the Traction threads. The Traction Threads come connected to the Sling from the manufacturer.
- 6.- Screwdriver: The Screwdriver is used to tighten the Fixation Screw (frontal screw of the Varitensor which fixes the Traction Threads).
- 7.- Passers: The detachable, single-use, retropubic Passers are used during the surgical procedure to drive the Traction Threads through to the abdominal wall.

**CONTRAINDICATIONS**

Do not implant the device:

- in patients with pre-existing infections or others condition that pose an unacceptable surgical risk as determined by the Physician, including comorbidities that may compromise tissue healing.
- in patients undergoing anticoagulation treatments.
- in patients who suffer urinary tract infections or obstructions.
- in pregnant women or women who are considering future pregnancies.
- in patients with renal insufficiency.
- in patients with immune diseases.
- in patients who are still growing or less than 18 years of age.
- in sub-urethral or peri-urethral areas with signs of tissue necrosis.
- in patients with known sensitivity or allergy to device materials.
- in patients undergoing concomitant bowel surgery.

WARNINGS AND PRECAUTIONS

Pre-operative:

- The patient should be counseled that there are alternatives to incontinence mesh surgeries, and the reason for choosing an adjustable mesh sling should be explained. Such alternatives include but are not limited to: Biofeedback, pelvic floor exercises (Kegel exercises), behavioral modification, pessary/vaginal tampons, transurethral bulking agents, medication, electrical stimulation and/or pads. The patient may also be offered a non-adjustable sling or alternatively, non-mesh surgical options (e.g.: fascial sling, stitches).
- Conservative, non-invasive treatments should be recommended before choosing a suburethral sling implant.
- Obtain patient consent of surgery with an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the device to be implanted is a permanent implant, and that some complications associated with the implant may require additional surgery that may or may not correct the complication. A multi-disciplinary team of physicians is recommended for patient counseling.
- The patient should be counseled that future pregnancies may negate the effects of surgical procedure and the patient may again become incontinent.
- The physician must determine the risks that should be included in the patient informed consent.
- Users should be familiar with the surgical technique for bladder neck suspensions and should be adequately trained in the placement Remeex System implantation procedure before employing the device.
- Surgical risks not associated directly with the use of the device, its function, or method of implantation require understanding and consideration by a qualified surgeon who is knowledgeable of anatomy, medical treatments and procedures for conditions appropriately treated with this device.
- The user must be qualified to manage any complications or adverse events associated to the device and the surgical procedure.
- As the anatomy of individual patients may vary greatly, for each procedure it is important that the intended plans for device advancement and the intended location for mesh placement are planned and known for each individual patient. Employment of imaging methods before and after mesh placement may aid in proper mesh placement and confirm absence of injury to non-target anatomical structures.
- The risks and benefits of performing the Remeex System procedure in patients with the following conditions should be carefully considered by the physician:
 - Patients with other pelvic floor meshes;
 - Overweight patients;
 - Concomitant pelvic organ surgery, including additional SUI treatment(s);
 - Patients with conditions that would compromise healing (in patients who suffer cancers of the vagina, cervix, or uterus or in patients who have received radiation in the area of treatment, etc.), and;
 - Patients with mild to light SUI.
- Remeex System device is not recommended for women in the first two years after their last pregnancy.
- Special care should be taken in patients with conditions that may be aggravated by their placement in modified dorsal lithotomy position.
- Special care should be taken in patients with anterior prolapse because of anatomical distortion of the pelvic floor.
- This product is provided sterile and must be sterile before its use. Before opening the product, the package must be inspected to ensure it is not damaged and that the sterility is not compromised.
- Do not reuse or resterilize this device. No component of the Female Remeex System is reusable. Reuse of the device could lead to potentially serious consequences for the patient that may lead to medical intervention, permanent injury or death. Reject and do not use any unpackaged system.
- The date of the last revision of this document is included in these Instructions For Use. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Intra-operative:

- This product is an implant that must be implanted in an operating theatre. Use strict aseptic and antibiotic coverage measures during the surgical procedure. Prophylactic antibiotics can be administered according to the physician's usual practice.
- Good surgical practices should be followed for management of contaminated or infected wounds with particular attention to the possibility of product infection.
- Cystoscopy should be performed to confirm bladder integrity after the Passer needle passage, as well as a digital rectal exam for detection of rectal perforation.
- Avoid excess tension of the mesh to minimize urinary retention/obstruction, bowel obstruction or the possibility of damage to tissues or internal organs. Mesh must not be placed in direct contact with bowel or visceral organs, including urinary bladder.
- Intra-operative bleeding may occur.
- Take special care during surgical intervention to avoid perforation of bladder, bowel, urethra, vaginal wall, intestine and other internal pelvic structures as well as damage of nerves, blood vessels, pubic bone, or mechanical disruption of the tissue. Attention to local anatomy and proper passage of needles will minimize risks.
- Do not contact any part of the Sling mesh with any staples, clips, or clamps as mechanical damage to the mesh may occur.
- No data is available on the use of the Sling with tissue adhesives and is not recommended.
- Omitting the vaginal and skin incisions to insert the device and Passers may result in excessive force and could cause injury to the patient.
- Surgeon(s) should change their gloves at each surgical site (gloves worn for the vaginal incision and sling placement should be changed when moving to the suprapubic surgical site).
- Do not turn the Manipulator without complete control. It is crucial to record the number of turns, the direction of the turns (both clockwise or counter-clockwise) and the date that it has been performed.
- While the external Manipulator is connected to the Varitensor, take the appropriate aseptic techniques and conditions to avoid contamination and infection at the site. The appropriate cleaning of the abdominal site must be performed by the appropriate medical personnel. During this time when the Manipulator and Varitensor are connected, the abdominal area must be protected using aseptic technique while covering the surgical field with appropriate sterile materials.
- To disconnect the external Manipulator from the Varitensor, gently pull on the Manipulator while rotating the Disconnecter until the point mark of the Disconnecter is aligned with the line mark of the Manipulator. Do not apply excessive pressure during this procedure.
- Clamps should only be used on the Traction Thread ends. The use of clamps could damage or weaken the strength of the Traction Threads. Damage to the Traction Threads between the Sling and the ends could cause system failure.
- The Passers, Manipulator, Disconnectors and Screwdriver are not permanent implants.
- During the readjustment procedure, take special care to avoid damaging the Traction Threads. Damage to the Traction Threads could cause system failure.
- The use of electrocautery is not recommended during the readjustment procedure given that the device material is sensitive to excessive temperatures and may cause further risks to the patient, product and/or healthcare specialist.
- The loosening and falling of the Fixation Screw into the surgical field could occur if surgeon mistakenly unscrews instead of tightening. Without the Fixation Screw, the Traction Threads cannot be fixed to the Varitensor.
- The Manipulator must not be connected to the Varitensor more than 3 consecutive days.
- The discarded device components (excess Traction Thread, Manipulator, Disconnectors, Screwdriver and Passers) may be a potential biohazard. The Passers have sharp tips, which could potentially protrude when not in use. Handle and dispose of the Passers, as well as the other discarded components in accordance with any local and federal laws regarding medical waste and disposal requirements for sharp devices, to prevent unnecessary injury or harm.

Post-operative:

- Instruct the patient to report dysuria, hematuria, paresthesia, cystoliths or any other abnormal symptoms to the surgeon as soon as possible.
- Post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient. The patient must be instructed to contact the surgeon immediately if any hemorrhaging occurs.
- Patient should report incontinence, voiding dysfunction or problems defecating to the surgeon as soon as possible.
- Traction Threads may fracture in case of overload due to overexertion of the patient.
- It is recommended that patients do not overly exert themselves physically (e.g.: heavy lifting) or partake in any moderate to high intensity sporting activity (e.g.: biking, jogging, etc.) for at least a minimum of two months. It is also recommended that patients avoid sexual intercourse during the first two months after surgery.
- The device is provided with multiple, product ("hospital") labels of which must be included in the patient medical file as well as an implant card to be administered to the patient for easy identification of the implant.
- Taking into account the condition of each individual patient, the surgeon will plan for post-operative follow-up. Follow-up is recommended at short/mid-term and long-term if possible.
- Sling tension-level readjustments may be made throughout the lifetime of the patient. However, wound dehiscence and/or infection may occur if the surgeon performs an incision for readjustment before the previous incision wound has healed completely.
- The Physician(s) and patient should report any serious incident in relation to Remeex System for Female Urinary Incontinence to the manufacturer and the national health authority having jurisdiction in their locale.

POTENTIAL ADVERSE REACTIONS

Known risks of surgical procedures for the treatment of urinary incontinence with a mid-urethral sling include:

- Temporary surrounding tissue irritation, allergic reaction or adverse reactions to foreign bodies may occur.
- Inflammation/swelling, adhesion formation, fistulas, abscess, seromas and/or hardening of the mesh caused by fibrosis.
- Bleeding (hemorrhage, perforation of vessels, hematuria, edema, erythema, hematoma).
- Perforation of bladder, bowel, vaginal wall, urethra, ureter, intestine and other pelvic structures.
- As a foreign body, the device could stimulate a pre-existing infection or suffer rejection.
- Extrusion, exposure, and/or migration of the mesh from the desired location. The device may also produce erosion through the vagina, urethra, bladder or other surrounding tissues or structures. Mesh contraction (shrinkage) is associated with vaginal tightening, vaginal shortening and vaginal pain.
- Urinary problems such as voiding symptoms/de novo incontinence in the form of urgency/overactive bladder (OAB), frequency or nocturia may appear after the intervention.
- Failure of the procedure resulting in recurrence (continued or worsening urinary incontinence).
- Altered bladder and/or bowel function such as urinary retention/obstruction or problems defecating.
- Infection, urinary infection.
- Acute and/or chronic pain, dysuria or other pains. Discomfort.
- Nerve damage, neuromuscular problems.
- Vaginal spotting, bleeding, abnormal discharge or scarring.
- Dyspareunia. Sexual partners may experience discomfort, irritation and/or pain during sexual activity.
- Foreign body sensitivity or palpability.

As with all surgical interventions, the patient is at risk of:

- Wound site complications (i.e.: excess granulation formation, wound dehiscence).
- Scarring.
- Physical deformity.
- Emotional problems.

In cases in which the Varitensor and/or the mesh need to be removed in part or whole, significant dissection may be required.

The occurrence of adverse reactions may require revisions and/or realization of specific, additional treatments, including but not limited to: antibiotics, blood transfusions, intravenous (IV) therapy, drainage, multiple restoration surgeries, removal of the device.

SURGICAL TECHNIQUE

Before using Remeex System for Female Urinary Incontinence the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical procedure to implant the Remeex System for Female Urinary Incontinence can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

Surgical technique:

Preparation of the implanted site:

1.- Abdominal Incision: With the patient in the lithotomic position, make a small transversal incision at the abdominal midline, 2 cm over the pubis. Then dissect the fat tissue until the fascia of the abdominal rectal muscles is clearly visible.

2.- Vaginal Incision: Perform a longitudinal incision at the midline of the anterior vaginal wall at the mid-urethral level. Dissect the paraurethral spaces to prepare the passage of the Passer needle to the abdominal site.

3.- Insertion of the Remeex Passer needles through the Retzius space: Remeex Passers are necessary for device delivery. They are used to drive the Traction Threads through the fascia of the abdominal rectal muscles. Insert the Remeex Passer needles in the para-urethral spaces, through the Retzius space, maintaining the tips away from the urethra and in continuous contact with the posterior wall of the pubis, until the tip of the Remeex Passer can be clearly seen at the abdominal incision. The Remeex Passers may be connected to the detachable handle which is used to help drive the needles and succeed in crossing the fascia of the abdominal rectal muscles. Bladder perforation may be avoided by leaving 3cm of distance between both needle tips when crossing the fascia of the abdominal rectal muscles. Take precautions to avoid going too lateral or too posterior in order to avoid blood vessels, bowel or other structural damage. Bladder mobilization techniques may also be used to avoid bladder perforation.

4.- Cystoscopy: Perform a cystoscopy to check the bladder integrity, and assure that the bladder is not perforated. If perforated, repair as necessary, withdraw the Remeex Passer and repeat step 3.

Placement of the Remeex System for Female Urinary Incontinence:

1.- Placement of the Sling and Traction Threads: Position the Sling at the midline of the anterior vaginal wall incision, at the level of the mid urethra. Then place the Traction Threads of the Sling through the holes in the proximal end of each Remeex Passer needle. Move the Remeex Passers forward by pushing them through the Retzius space towards the abdominal site

incisions or pulling them from the abdominal site until the Traction Threads appear at the abdominal incision. The surgical assistant should clamp the Traction Thread ends with a mosquito clamp. Do not clamp the Traction Threads in any other area apart from the tip. Fix the Sling with two absorbable suture knots to avoid Remeex Sling antero-posterior or lateral displacement of the Sling. The Sling must be in full contact with the suburethral tissue.

2.- Varitensor placement: Before starting this part of the surgery, the surgeon and the assistant need to change the surgical gloves to avoid contaminating the Varitensor with the vaginal flora. Insert the Traction Threads through the Varitensor's respective lateral reception holes, so that the Traction Threads come out through the central outcome hole (at the Varitensor midline). Before inserting the Traction Threads through the Varitensor, align the marks of the external ring of the Fixation Screw by rotating the Manipulator (the marks need to be horizontally aligned with the Traction Thread reception holes). Once this step is completed, have the surgical assistant hold the Varitensor prosthesis in the midline, 4 cm over the fascia of the abdominal rectal muscles, and then fix the Traction Threads by tightening the Fixation Screw. Verify that the same length of Traction Thread has been left at each side between the Varitensor and the fascia. Cut and discard the excess Traction Thread protruding from the Fixation Screw. Rotate the Manipulator clockwise, winding the Traction Threads into the Varitensor, which will lower the Varitensor until they rest over the fascia of the abdominal rectus muscles, where it will remain for future adjustment. During the Manipulator rotation operation, it is very important to maintain a slight upward pull on the Varitensor, using the other hand. This minimal tension is optimal for permitting the Traction Threads to be wound properly into the Varitensor.

3.- Vaginal and abdominal incision closure: Close the abdominal incision with the Manipulator protruding perpendicular to the abdominal wall. Close the vaginal incision.

Post surgical regulation and Manipulator disconnection:


The following morning, fill the bladder with 300 c.c. of saline through a urethral catheter and ask the patient to stand up and perform Valsalva maneuvers. If necessary, rotate the Manipulator clockwise checking continence level every 4 complete turns, until continence is achieved. Invite the patient to urinate and measure the residual in the bladder after urination. If residual is less than 100 c.c. disconnect the Manipulator and discharge the patient. If residual is more than 100 c.c., decrease the Sling tension by rotating the Manipulator counterclockwise, helping to descend the urethra using a rigid probe. To disconnect the Manipulator from the Varitensor, insert the Disconnecter inside the Manipulator and rotate the Disconnecter $\frac{1}{4}$ of a turn in relation to the Manipulator while slightly pulling the Manipulator. *It is important to explain to the patient that even though she is leaving the hospital with continence, some degree of incontinence may re-appear during the first month after the intervention. If this is the case, it is possible to re adjust the Sling level by means of a minor surgical procedure.*

Medium-Long term ability to readjust the sling support level:

(For medium-long term readjustment, wait at least 30 days from hospital discharge to ensure that the patient is completely recovered from surgery).

Remeex System for Female Urinary Incontinence offers the possibility to modify the Sling support level whenever needed throughout the patient's lifetime. During the month following the surgical procedure the tissue inflammation will diminish considerably, the manipulated organs will return to their normal anatomical position and the patient will return to her normal physical activity level. Some patients may show some degree of incontinence, in which case the Sling can be readjusted to restore continence by performing a minor surgical procedure. This is accomplished using the single-use, pre-sterilized Remeex Readjustment Kit which contains 1 Manipulator and 2 Disconnectors.

Sling Level Readjustment: In an appropriate surgical environment, after injecting local anaesthesia, make a minimal incision at the abdominal midline, just over the pubis in order to locate the Varitensor. Connect the sterile Manipulator to the implanted Varitensor and close the incision with the Manipulator protruding from the abdomen. The Sling readjustment is performed with the patient in standing position, performing the normal pressure maneuvers that cause the patient to become incontinent. Fill the bladder with 300 c.c. of saline through a urethral catheter and ask the patient to stand up and perform Valsalva maneuvers. If necessary, rotate the Manipulator clockwise checking continence level every 4 complete turns, until continence is achieved (there is no leakage). Invite the patient to urinate and measure the residual in the bladder after urination. If residual is less than 100 c.c. disconnect the Manipulator from the Varitensor, close the incision and discharge the patient. If residual is more than 100 c.c., decrease the Sling tension by rotating the Manipulator counterclockwise, helping to descend the urethra using a rigid probe. To disconnect the Manipulator from the Varitensor, insert the Disconnecter inside the Manipulator and rotate the Disconnecter $\frac{1}{4}$ of a turn in relation to the Manipulator while slightly pulling the Manipulator. Close the abdominal incision.


Remeex System for Female Urinary Incontinence removal:

The removal technique of Remeex System for Female Urinary Incontinence can be performed under regional or general anesthesia. It is advisable to use the hospital antibiotic protocol for implants.

- 1.- Abdominal incision: With the patient in modified dorsal lithotomy position, prepare the surgical field. Start by making a 4cm transversal abdominal incision 2 fingers over the symphysis pubis. Identify the Varitensor by palpation and/or with a long anesthetic needle. Dissect tissue until the Varitensor is fully visible. Connect the Manipulator to the Varitensor and unwind the internal traction thread, cut the thread and discard the Varitensor. Rinse the entire surgical site with antibiotic saline solution and close the abdominal incision.
- 2.- Vaginal incision: With the patient in modified dorsal lithotomy position, prepare the surgical field with a vaginal valve. Start by making a longitudinal incision in the anterior vaginal wall, preparing the submucosal space to present the best possible view of the urethral tract.
- 3.- Identify the mesh by direct view and palpation. Dissect the mesh from the surrounding structures. Avoid cutting the mesh whenever possible. Dissection of the mesh should continue by slightly pulling it to identify the mesh path. When mesh is completely dissected, remove the mesh and the traction threads from the vaginal side. Rinse the entire surgical site with antibiotic saline solution and close the vaginal incision.

MECHANICAL FUNCTION STEPS

- Figure 1.- Horizontal alignment of the external ring marks of the Fixation Screw by rotating the Manipulator.
 Figure 2.- Passing of the Traction Threads through the Varitensor.
 Figure 3.- Traction Threads fixation with Fixation Screw.
 Figure 4.- Cutting of the excess Traction Threads.
 Figure 5.- Sling level regulation by rotating the Manipulator.
 Figure 6.- Disconnecter insertion into the Manipulator.
 Figure 7.- 90° rotation of the Disconnecter in relation to the Manipulator.
 Figure 8.- Disconnection of the Manipulator.

MECHANICAL CHARACTERISTICS /SPECIFICATIONS

External manipulator function: Every 360° turn of the external Remeex Manipulator causes an increase or decrease of the Sling's suspension level by approximately 1mm.

Wind up capacity: The wind up capacity of the Varitensor is 15 cm of Traction Threads per each side.

Maximum traction stress: The Varitensor can support the force of 18 Newtons produced by the Traction Threads, with a minimum safety coefficient of 5.

MATERIALS USED IN MANUFACTURING PROCESSPatient contacting materials:

Varitensor and Manipulator: Ultra-high molecular weight Polyethylene (UHMWPE), Titanium

Sling: Polypropylene monofilament mesh (colorless), blue polyvinylidene fluoride (PVDF) monofilament Traction Threads (colorant: solvent blue 104, CAS no.: 116-75-6)

Passers*: Stainless steel

No phthalates, latex or other hazardous materials are used in the manufacture of the any Remeex System implant components.

*Stainless steel components contain nickel

STERILIZATION

The device is sterilized by Ethylene Oxide. This product is a single use device and it must not be re-sterilized or re-used. Do not use if the packaging has been damaged.

STORAGE

It is recommended to store the product in a dry place at room temperature, away from heat and sunlight.

Do not use after the expiration date detailed on the package.

SINGLE USE PRODUCT DISCLAIMER

Device intended to be used only once and for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

SYMBOLS USED ON LABELING

CATALOGUE NUMBER

DISTRIBUTOR

STERILIZED USING ETHYLENE OXIDE

SERIAL NUMBER

CONSULT INSTRUCTIONS FOR USE

SINGLE STERILE BARRIER SYSTEM

UNIQUE DEVICE IDENTIFIER

USE BY DATE

DATE OF IMPLANTATION

MEDICAL DEVICE

MANUFACTURER

PATIENT NAME OR PATIENT ID

DO NOT REUSE

KEEP DRY

NAME/ADDRESS OF IMPLANTING HEALTHCARE INSTITUTION/PROVIDER

DO NOT RESTERILIZE

DATE OF MANUFACTURE

INFORMATION WEBSITE FOR PATIENTS

DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE

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